

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | rp98927       | Batch NO. :          | MVM2500713              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 56 Min 07 Sec   |
| Total Counter :  | 550226        | Total Running Time : | 862 Hour 3 Min          |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 52                | %       | Manual ;           | 11           | RPM        |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm    |
| Target Quantity :                  |             | 2000000           | Tablets |                    | Max :        | 12.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 30                | Min     | Supply Time ;      | 0.2          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 2.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 15.52             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 14.92             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 12.18             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 11.94             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 11.70             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 8.95              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 8.36              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.15       | kN  | Min : | 2.12       | kN  | Punch NO.; | 43 |
| SD :                 | 0.90       | %   | Max : | 2.20       | kN  | Punch NO.; | 11 |
| Mean Main Pressure : | 12.00      | kN  | Min : | 11.00      | kN  | Punch NO.; | 40 |
| SD :                 | 3.60       | %   | Max : | 13.20      | kN  | Punch NO.; | 15 |
| Filling Depth :      | Increase ; | 585 | Times | Decrease ; | 546 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 11       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.67    | mm  | Main Thickness ; | 1.68     | mm           |
| Filling Depth :               | 9.08    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.802   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 14      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.25 | kN |             |    |
|                  | Max ;     | 0.32 | kN | Max Punch ; | 17 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 46 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 550226             | Tablets |       |        |         |
| A.W.C Counter : | 548835             | Tablets |       |        |         |
| HEP :           | 77 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 548758<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | rp98927       | Batch NO. :          | MVM2500713              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 56 Min 07 Sec   |
| Total Counter :  | 550226        | Total Running Time : | 862 Hour 3 Min          |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 55                | %       | Manual ;           | 11           | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |        | 2000000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 15.37             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 14.78             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 12.06             | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 11.82             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 11.59             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 8.87              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 8.28              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.08       | kN  | Min : | 2.05       | kN  | Punch NO.; | 2  |
| SD :                 | 0.70       | %   | Max : | 2.11       | kN  | Punch NO.; | 8  |
| Mean Main Pressure : | 11.94      | kN  | Min : | 10.69      | kN  | Punch NO.; | 37 |
| SD :                 | 5.40       | %   | Max : | 13.33      | kN  | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 226 | Times | Decrease ; | 244 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.68    | mm  | Main Thickness ; | 1.69     | mm           |
| Filling Depth :               | 8.91    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.802   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 14      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.25 | kN |             |    |
|                  | Max ;     | 0.32 | kN | Max Punch ; | 17 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 46 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 550226             | Tablets |       |        |         |
| A.W.C Counter : | 548835             | Tablets |       |        |         |
| HEP :           | 77 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 548758<br>(100.0%) | Tablets |       |        |         |