

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | pb97604       | Batch NO. :          | MVM2500713              |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 52 Min 08 Sec   |
| Total Counter :  | 1100494       | Total Running Time : | 865 Hour 59 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 52                | %       | Manual ;           | 11    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 2000000           | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 15.90             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 15.29             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 12.47             | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 12.23             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 11.98             | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 9.17              | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 8.56              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.18       | kN  | Min : | 2.14       | kN  | Punch NO.; | 17 |
| SD :                 | 0.80       | %   | Max : | 2.22       | kN  | Punch NO.; | 27 |
| Mean Main Pressure : | 12.21      | kN  | Min : | 11.18      | kN  | Punch NO.; | 40 |
| SD :                 | 3.80       | %   | Max : | 13.24      | kN  | Punch NO.; | 30 |
| Filling Depth :      | Increase ; | 660 | Times | Decrease ; | 607 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 11       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.67   | mm  | Main Thickness ; | 1.68     | mm           |
| Filling Depth :               | 9.17   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.786  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 22     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.35 | kN | Max Punch ; | 28 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 46 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1100494             | Tablets |       |            |         |
| A.W.C Counter : | 1097535             | Tablets |       |            |         |
| HEP :           | 190 (0.0%)          | Tablets | LEP : | 237 (0.0%) | Tablets |
| GOOD :          | 1097108<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MEBEV 135<br>mgSTMV2A00 |
| User ID :        | pb97604       | Batch NO. :          | MVM2500713              |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 52 Min 08 Sec   |
| Total Counter :  | 1100494       | Total Running Time : | 865 Hour 59 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 55                | %       | Manual ;           | 11    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 2000000           | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30                | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 16.19             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 15.57             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 12.70             | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 12.45             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 12.20             | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 9.34              | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 8.72              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.27       | kN  | Min : | 2.23       | kN  | Punch NO.; | 2  |
| SD :                 | 0.80       | %   | Max : | 2.33       | kN  | Punch NO.; | 27 |
| Mean Main Pressure : | 12.49      | kN  | Min : | 11.21      | kN  | Punch NO.; | 19 |
| SD :                 | 4.80       | %   | Max : | 13.65      | kN  | Punch NO.; | 28 |
| Filling Depth :      | Increase ; | 424 | Times | Decrease ; | 443 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.68   | mm  | Main Thickness ; | 1.69     | mm           |
| Filling Depth :               | 9.06   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.786  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 22     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 14     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.27 | kN |             |    |
|                  | Max ;     | 0.34 | kN | Max Punch ; | 17 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 46 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1100494             | Tablets |       |            |         |
| A.W.C Counter : | 1097535             | Tablets |       |            |         |
| HEP :           | 190 (0.0%)          | Tablets | LEP : | 237 (0.0%) | Tablets |
| GOOD :          | 1097108<br>(100.0%) | Tablets |       |            |         |