

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Prop 40mg STPJ2B00    |
| User ID :        | pb97604       | Batch NO. :          | PJM2500240            |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 01 Min 25 Sec |
| Total Counter :  | 1750980       | Total Running Time : | 897 Hour 25 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 35       | %       | Manual ;           | 11           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 31.14    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 29.94    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 24.43    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 23.95    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 23.47    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 17.96    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 16.77    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.84         | kN    | Min :      | 1.78  | kN    | Punch NO.; | 18 |
| SD :                 | 1.50         | %     | Max :      | 1.91  | kN    | Punch NO.; | 46 |
| Mean Main Pressure : | 23.89        | kN    | Min :      | 22.73 | kN    | Punch NO.; | 26 |
| SD :                 | 2.30         | %     | Max :      | 24.92 | kN    | Punch NO.; | 2  |
| Filling Depth :      | Increase ; 7 | Times | Decrease ; | 49    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 159000            | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 9                 | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.55   | mm  | Main Thickness ; | 1.56              | mm           |
| Filling Depth :               | 5.49   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.818  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 15     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 8      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.35 | kN |             |    |
|                  | Max ;     | 0.42 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 34 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1750980             | Tablets |       |        |         |
| A.W.C Counter : | 1746826             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1746826<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Prop 40mg STPJ2B00    |
| User ID :        | pb97604       | Batch NO. :          | PJM2500240            |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 01 Min 25 Sec |
| Total Counter :  | 1750980       | Total Running Time : | 897 Hour 25 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 35        | %       | Manual ;           | 11    | RPM                  |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min : | 2.50 mm              |
|                                    |        | 4.5-10mm  |         |                    | Max : | 8.00 mm              |
| Target Quantity :                  |        | 2000000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.600     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20        | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 29.90     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 28.75     | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 23.46     | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 23.00     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 22.54     | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 17.25     | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 16.10     | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |           |         |                    |       |                      |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.98  | kN    | Min :      | 1.96  | kN    | Punch NO.; | 24 |
|                      |            | 0.50  | %     | Max :      | 2.01  | kN    | Punch NO.; | 38 |
| Mean Main Pressure : |            | 23.04 | kN    | Min :      | 21.01 | kN    | Punch NO.; | 40 |
|                      | SD :       | 4.00  | %     | Max :      | 25.10 | kN    | Punch NO.; | 34 |
| Filling Depth :      | Increase ; | 23    | Times | Decrease ; | 33    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.55   | mm  | Main Thickness ; | 1.56     | mm           |
| Filling Depth :               | 5.51   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.818  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 15     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.35 | kN |             |    |
|                  | Max ;     | 0.42 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 22 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1750980             | Tablets |       |        |         |
| A.W.C Counter : | 1746826             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1746826<br>(100.0%) | Tablets |       |        |         |