

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METF 500MG STEW1A00
User ID :	pb97604	Batch NO. :	EWM2505650
Print Interval :	50000 Tablets	Running Time :	06 Hour 45 Min 28 Sec
Total Counter :	1200440	Total Running Time :	959 Hour 17 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	107	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		4500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		180	Min	Supply Time ;	0.8	Sec
Lower Punch Interval Time(S2) ;		150	Min	Supply Time ;	2.2	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	15.75	kN	No		
HEP :	25.0 %	15.14	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	12.35	kN	Yes 10	Times	
Ref :		12.11	kN			
LCP :	2.0 %	11.87	kN	Yes 10	Times	
LEP :	25.0 %	9.08	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	8.48	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.86	kN	Min :	2.81	kN	Punch NO.;	16
SD :	0.60	%	Max :	2.91	kN	Punch NO.;	24
Mean Main Pressure :	11.99	kN	Min :	10.25	kN	Punch NO.;	36
SD :	5.20	%	Max :	13.36	kN	Punch NO.;	6
Filling Depth :	Increase ;	327	Times	Decrease ;	305	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	28	RPM	Capacity :	178080	Tablets/hour
Feeder : Status ;	Auto		Speed ;	30	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.35	mm	Pre Thickness ;	3.36	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	2.41	mm
Filling Depth :	7.74	mm	Current Cam :	Cam B	
Main Air pressure :	0.803	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	132	Min			
Lower Punch Remain Time (S2):	113	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.01	kN	Max Punch ;	36
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.13	kN	Max Punch ;	51
Ejection Force :	Average ;	0.08	kN		
	Max ;	0.15	kN	Max Punch ;	46

Tablet Data (N/A)

Total Counter :	1200440	Tablets			
A.W.C Counter :	1192792	Tablets			
HEP :	245 (0.0%)	Tablets	LEP :	36 (0.0%)	Tablets
GOOD :	1192511 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METF 500MG STEW1A00
User ID :	pb97604	Batch NO. :	EWM2505650
Print Interval :	50000 Tablets	Running Time :	06 Hour 45 Min 28 Sec
Total Counter :	1200440	Total Running Time :	959 Hour 17 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	100	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		4500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		180	Min	Supply Time ;	0.8	Sec
Lower Punch Interval Time(S2) ;		150	Min	Supply Time ;	2.2	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	15.84	kN	No		
HEP :	25.0 %	15.23	kN	Yes 10	Rotations 10	Tablets
HCP :	2.0 %	12.43	kN	Yes 10	Times	
Ref :		12.18	kN			
LCP :	2.0 %	11.94	kN	Yes 10	Times	
LEP :	25.0 %	9.14	kN	Yes 10	Rotations 10	Tablets
LSP :	30.0 %	8.53	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.80	kN	Min :	2.73	kN	Punch NO.;	30
SD :	0.70	%	Max :	2.85	kN	Punch NO.;	9
Mean Main Pressure :	12.26	kN	Min :	10.24	kN	Punch NO.;	36
SD :	6.60	%	Max :	14.16	kN	Punch NO.;	6
Filling Depth :	Increase ;	277	Times	Decrease ;	283	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	28	RPM	Capacity :	178080	Tablets/hour
Feeder : Status ;	Auto		Speed ;	28	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	3.31	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	5.40	mm	Main Thickness ;	2.41	mm
Filling Depth :	7.44	mm	Current Cam :	Cam B	
Main Air pressure :	0.803	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	132	Min			
Lower Punch Remain Time (S2):	113	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.02	kN	Max Punch ;	36
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.13	kN	Max Punch ;	51
Ejection Force :	Average ;	0.08	kN		
	Max ;	0.15	kN	Max Punch ;	46

Tablet Data (N/A)

Total Counter :	1200440	Tablets			
A.W.C Counter :	1192792	Tablets			
HEP :	245 (0.0%)	Tablets	LEP :	36 (0.0%)	Tablets
GOOD :	1192511 (100.0%)	Tablets			