

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | BUME 1MG STJU1A00     |
| User ID :        | am98755       | Batch NO. :          | JUM2500301            |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 08 Min 33 Sec |
| Total Counter :  | 300134        | Total Running Time : | 982 Hour 17 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |           |
|------------------------------------|--------|-----------|---------|--------------------|-------|-----------|
| Feeder :                           | Auto ; | 13        | %       | Manual ;           | 15    | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min : | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max : | 8.00 mm   |
| Target Quantity :                  |        | 360000    | Tablets |                    |       |           |
| Air Pressure Low Limit :           |        | 0.600     | MPa     |                    |       |           |
| Oil Lubrication :                  |        |           |         |                    |       |           |
| Upper Punch Interval Time(S1) ;    |        | 15        | Min     | Supply Time ;      | 2.0   | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15        | Min     | Supply Time ;      | 2.0   | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |       |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |       |           |
|                                    |        |           |         | Stop Condition ;   |       |           |
| HSP :                              | 30.0 % | 18.28     | kN      | No                 |       |           |
| HEP :                              | 25.0 % | 17.58     | kN      | Yes                | 15    | Rotations |
| HCP :                              | 2.0 %  | 14.34     | kN      | Yes                | 20    | Times     |
| Ref :                              |        | 14.06     | kN      |                    |       |           |
| LCP :                              | 2.0 %  | 13.78     | kN      | Yes                | 20    | Times     |
| LEP :                              | 25.0 % | 10.55     | kN      | Yes                | 15    | Rotations |
| LSP :                              | 30.0 % | 9.84      | kN      | No                 |       |           |
| SD limit :                         | 40.0 % |           |         | Yes                |       |           |
| Pre HSP                            |        | 35.00     | kN      | Yes                |       |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations |
| Empty Punch No. :                  | None   |           |         |                    |       |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 4.54       | kN  | Min : | 4.23       | kN  | Punch NO.; | 42 |
| SD :                 | 3.00       | %   | Max : | 4.95       | kN  | Punch NO.; | 48 |
| Mean Main Pressure : | 11.82      | kN  | Min : | 10.48      | kN  | Punch NO.; | 42 |
| SD :                 | 5.00       | %   | Max : | 13.51      | kN  | Punch NO.; | 48 |
| Filling Depth :      | Increase ; | 329 | Times | Decrease ; | 295 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 140400   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 2        | RPM          |
| Pre Upper Punch Position ;    | 2.82    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 3.80    | mm  | Pre Thickness ;  | 0.98     | mm           |
| Main Upper Punch Position ;   | 3.49    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.01    | mm  | Main Thickness ; | 0.52     | mm           |
| Filling Depth :               | 4.58    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.798   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 12 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 14 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 27 |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 300134             | Tablets |       |           |         |
| A.W.C Counter : | 296231             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 40 (0.0%) | Tablets |
| GOOD :          | 296191<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | BUME 1MG STJU1A00     |
| User ID :        | am98755       | Batch NO. :          | JUM2500301            |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 08 Min 33 Sec |
| Total Counter :  | 300134        | Total Running Time : | 982 Hour 17 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |           |         |                    |       |           |
|------------------------------------|--------|-----------|---------|--------------------|-------|-----------|
| Feeder :                           | Auto ; | 13        | %       | Manual ;           | 15    | RPM       |
| Filling Cam :                      |        | Cam B     |         | Cam Range :        | Min : | 2.50 mm   |
|                                    |        | 4.5-10mm  |         |                    | Max : | 8.00 mm   |
| Target Quantity :                  |        | 360000    | Tablets |                    |       |           |
| Air Pressure Low Limit :           |        | 0.600     | MPa     |                    |       |           |
| Oil Lubrication :                  |        |           |         |                    |       |           |
| Upper Punch Interval Time(S1) ;    |        | 15        | Min     | Supply Time ;      | 2.0   | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15        | Min     | Supply Time ;      | 2.0   | Sec       |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |       |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |       |           |
|                                    |        |           |         | Stop Condition ;   |       |           |
| HSP :                              | 30.0 % | 19.78     | kN      | No                 |       |           |
| HEP :                              | 25.0 % | 19.02     | kN      | Yes                | 15    | Rotations |
| HCP :                              | 2.0 %  | 15.52     | kN      | Yes                | 20    | Times     |
| Ref :                              |        | 15.21     | kN      |                    |       |           |
| LCP :                              | 2.0 %  | 14.91     | kN      | Yes                | 20    | Times     |
| LEP :                              | 25.0 % | 11.41     | kN      | Yes                | 15    | Rotations |
| LSP :                              | 30.0 % | 10.65     | kN      | No                 |       |           |
| SD limit :                         | 40.0 % |           |         | Yes                |       |           |
| Pre HSP                            |        | 35.00     | kN      | Yes                |       |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |       |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 4.63       | kN  | Min : | 4.38       | kN  | Punch NO.; | 26 |
| SD :                 | 3.10       | %   | Max : | 4.95       | kN  | Punch NO.; | 17 |
| Mean Main Pressure : | 15.20      | kN  | Min : | 13.87      | kN  | Punch NO.; | 26 |
| SD :                 | 4.00       | %   | Max : | 16.74      | kN  | Punch NO.; | 17 |
| Filling Depth :      | Increase ; | 253 | Times | Decrease ; | 232 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 140400   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 2        | RPM          |
| Pre Upper Punch Position ;    | 2.80    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 3.80    | mm  | Pre Thickness ;  | 1.00     | mm           |
| Main Upper Punch Position ;   | 3.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.01    | mm  | Main Thickness ; | 0.51     | mm           |
| Filling Depth :               | 3.74    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.798   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 7       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 12 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 53 |
| Ejection Force : | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 27 |

### Tablet Data (N/A)

|                 |                    |         |       |           |         |
|-----------------|--------------------|---------|-------|-----------|---------|
| Total Counter : | 300134             | Tablets |       |           |         |
| A.W.C Counter : | 296231             | Tablets |       |           |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 40 (0.0%) | Tablets |
| GOOD :          | 296191<br>(100.0%) | Tablets |       |           |         |