

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 400MGSTIU2B00
User ID :	pb97604	Batch NO. :	IUM2500340
Print Interval :	50000 Tablets	Running Time :	00 Hour 51 Min 35 Sec
Total Counter :	100052	Total Running Time :	986 Hour 20 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	20.32	kN	No		
HEP :	25.0 %	19.54	kN	Yes	15	Rotations
HCP :	2.0 %	15.94	kN	Yes	20	Times
Ref :		15.63	kN			
LCP :	2.0 %	15.32	kN	Yes	20	Times
LEP :	25.0 %	11.72	kN	Yes	15	Rotations
LSP :	30.0 %	10.94	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	5.11	kN	Min :	4.80	kN	Punch NO.;	18
SD :	2.40	%	Max :	5.37	kN	Punch NO.;	8
Mean Main Pressure :	15.35	kN	Min :	14.32	kN	Punch NO.;	31
SD :	3.60	%	Max :	16.45	kN	Punch NO.;	14
Filling Depth :	Increase ;	39	Times	Decrease ;	68	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.45	mm	Main Thickness ;	1.46	mm
Filling Depth :	6.85	mm	Current Cam :	Cam C	
Main Air pressure :	0.803	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	13	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.10	kN		
	Max ;	0.23	kN	Max Punch ;	41
Lower Punch :	Average ;	0.13	kN		
	Max ;	0.17	kN	Max Punch ;	41
Ejection Force :	Average ;	0.06	kN		
	Max ;	0.08	kN	Max Punch ;	13

Tablet Data (N/A)

Total Counter :	100052	Tablets			
A.W.C Counter :	89601	Tablets			
HEP :	69 (0.1%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	89531 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBUPRO 400MGSTIU2B00
User ID :	pb97604	Batch NO. :	IUM2500340
Print Interval :	50000 Tablets	Running Time :	00 Hour 51 Min 35 Sec
Total Counter :	100052	Total Running Time :	986 Hour 20 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	200	%	Manual ;	15	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		2025000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		5.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	19.38	kN	No		
HEP :	25.0 %	18.63	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	15.21	kN	Yes 20	Times	
Ref :		14.91	kN			
LCP :	2.0 %	14.61	kN	Yes 20	Times	
LEP :	25.0 %	11.18	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	10.44	kN	No		
SD limit :	40.0 %			Yes		
Pre HSP		35.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	4.97	kN	Min :	4.76	kN	Punch NO.;	10
SD :	2.00	%	Max :	5.18	kN	Punch NO.;	50
Mean Main Pressure :	14.59	kN	Min :	13.37	kN	Punch NO.;	30
SD :	4.50	%	Max :	15.94	kN	Punch NO.;	53
Filling Depth :	Increase ;	56	Times	Decrease ;	48	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	50	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	2.51	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.60	mm	Main Thickness ;	1.61	mm
Filling Depth :	7.00	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.803	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	13	Min			
Lower Punch Remain Time (S2):	13	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.10	kN		
	Max ;	0.21	kN	Max Punch ;	43
Lower Punch :	Average ;	0.13	kN		
	Max ;	0.17	kN	Max Punch ;	41
Ejection Force :	Average ;	0.06	kN		
	Max ;	0.08	kN	Max Punch ;	13

Tablet Data (N/A)

Total Counter :	100052	Tablets			
A.W.C Counter :	89654	Tablets			
HEP :	69 (0.1%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	89584 (99.9%)	Tablets			