

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBUPRO<br>400MGSTIU2B00 |
| User ID :        | yd98913       | Batch NO. :          | IUM2500340              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 59 Min 42 Sec   |
| Total Counter :  | 1450678       | Total Running Time : | 995 Hour 29 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 200               | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 2025000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 20.95             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 20.15             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 16.44             | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 16.12             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 15.79             | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 12.09             | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 11.28             | kN      | No                 |              |           |
| SD limit :                         | 40.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 35.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 5.18       | kN  | Min : | 4.91       | kN  | Punch NO.; | 26 |
| SD :                 | 2.80       | %   | Max : | 5.56       | kN  | Punch NO.; | 47 |
| Mean Main Pressure : | 16.73      | kN  | Min : | 15.18      | kN  | Punch NO.; | 25 |
| SD :                 | 4.70       | %   | Max : | 18.36      | kN  | Punch NO.; | 8  |
| Filling Depth :      | Increase ; | 495 | Times | Decrease ; | 450 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 40       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.45    | mm  | Main Thickness ; | 1.46     | mm           |
| Filling Depth :               | 9.15    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.767   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 12      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 32 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 16 |
| Ejection Force : | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 1450678             | Tablets |       |          |         |
| A.W.C Counter : | 1435829             | Tablets |       |          |         |
| HEP :           | 706 (0.0%)          | Tablets | LEP : | 1 (0.0%) | Tablets |
| GOOD :          | 1435122<br>(100.0%) | Tablets |       |          |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBUPRO<br>400MGSTIU2B00 |
| User ID :        | yd98913       | Batch NO. :          | IUM2500340              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 59 Min 42 Sec   |
| Total Counter :  | 1450678       | Total Running Time : | 995 Hour 29 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |           |
|------------------------------------|--------|-------------------|---------|--------------------|-------|-----------|
| Feeder :                           | Auto ; | 200               | %       | Manual ;           | 15    | RPM       |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm   |
| Target Quantity :                  |        | 2025000           | Tablets |                    | Max : | 12.00 mm  |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |           |
| Oil Lubrication :                  |        |                   |         |                    |       |           |
| Upper Punch Interval Time(S1) ;    |        | 15                | Min     | Supply Time ;      | 2.0   | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 2.0   | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |           |
|                                    |        |                   |         | Stop Condition ;   |       |           |
| HSP :                              | 30.0 % | 21.73             | kN      | No                 |       |           |
| HEP :                              | 25.0 % | 20.89             | kN      | Yes                | 15    | Rotations |
| HCP :                              | 2.0 %  | 17.05             | kN      | Yes                | 20    | Times     |
| Ref :                              |        | 16.71             | kN      |                    |       |           |
| LCP :                              | 2.0 %  | 16.38             | kN      | Yes                | 20    | Times     |
| LEP :                              | 25.0 % | 12.53             | kN      | Yes                | 15    | Rotations |
| LSP :                              | 30.0 % | 11.70             | kN      | No                 |       |           |
| SD limit :                         | 40.0 % |                   |         | Yes                |       |           |
| Pre HSP                            |        | 35.00             | kN      | Yes                |       |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |       |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 5.14       | kN  | Min : | 4.73       | kN  | Punch NO.; | 33 |
| SD :                 | 4.50       | %   | Max : | 6.04       | kN  | Punch NO.; | 39 |
| Mean Main Pressure : | 16.98      | kN  | Min : | 14.94      | kN  | Punch NO.; | 35 |
| SD :                 | 5.60       | %   | Max : | 18.81      | kN  | Punch NO.; | 52 |
| Filling Depth :      | Increase ; | 390 | Times | Decrease ; | 333 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 40                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50    | mm  | Pre Thickness ;  | 2.51              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 4.60    | mm  | Main Thickness ; | 1.61              | mm           |
| Filling Depth :               | 8.39    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.767   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 12      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 16 |
| Ejection Force : | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                     |         |       |          |         |
|-----------------|---------------------|---------|-------|----------|---------|
| Total Counter : | 1450678             | Tablets |       |          |         |
| A.W.C Counter : | 1435829             | Tablets |       |          |         |
| HEP :           | 706 (0.0%)          | Tablets | LEP : | 1 (0.0%) | Tablets |
| GOOD :          | 1435122<br>(100.0%) | Tablets |       |          |         |