

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBUPRO<br>600MGSTIU2C00 |
| User ID :        | pm98770       | Batch NO. :          | IUM2501260              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 54 Min 11 Sec   |
| Total Counter :  | 205614        | Total Running Time : | 1001 Hour 51 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 365               | %       | Manual ;           | 19           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 1350000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 24.30             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 23.36             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 19.06             | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 18.69             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 18.32             | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 14.02             | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 13.08             | kN      | No                 |              |           |
| SD limit :                         | 40.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 35.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 7.74          | kN    | Min :         | 6.97  | kN | Punch NO.; | 14 |
| SD :                 | 5.40          | %     | Max :         | 8.55  | kN | Punch NO.; | 17 |
| Mean Main Pressure : | 18.44         | kN    | Min :         | 15.78 | kN | Punch NO.; | 43 |
| SD :                 | 5.50          | %     | Max :         | 20.34 | kN | Punch NO.; | 14 |
| Filling Depth :      | Increase ; 50 | Times | Decrease ; 62 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 66       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 4.30     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.35    | mm  | Main Thickness ; | 3.35     | mm           |
| Filling Depth :               | 11.34   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.773   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 8       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 1  |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 205614             | Tablets |       |        |         |
| A.W.C Counter : | 197464             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 197464<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | IBUPRO<br>600MGSTIU2C00 |
| User ID :        | pm98770       | Batch NO. :          | IUM2501260              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 54 Min 11 Sec   |
| Total Counter :  | 205614        | Total Running Time : | 1001 Hour 51 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 277               | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 1350000           | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 15                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 27.00             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 25.96             | kN      | Yes                | 15    | Rotations 20 Tablets |
| HCP :                              | 2.0 %  | 21.18             | kN      | Yes                | 20    | Times                |
| Ref :                              |        | 20.77             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 20.35             | kN      | Yes                | 20    | Times                |
| LEP :                              | 25.0 % | 15.58             | kN      | Yes                | 15    | Rotations 20 Tablets |
| LSP :                              | 30.0 % | 14.54             | kN      | No                 |       |                      |
| SD limit :                         | 40.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 35.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 7.98       | kN | Min : | 6.97       | kN | Punch NO.; | 30 |
| SD :                 | 5.50       | %  | Max : | 8.85       | kN | Punch NO.; | 12 |
| Mean Main Pressure : | 20.33      | kN | Min : | 17.73      | kN | Punch NO.; | 44 |
| SD :                 | 5.70       | %  | Max : | 22.96      | kN | Punch NO.; | 51 |
| Filling Depth :      | Increase ; | 36 | Times | Decrease ; | 10 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 18      | RPM | Capacity :       | 114480   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 50       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 4.30     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.45    | mm  | Main Thickness ; | 3.45     | mm           |
| Filling Depth :               | 11.35   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.773   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 8       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 52 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 41 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 16 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 205614             | Tablets |       |        |         |
| A.W.C Counter : | 197517             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 197517<br>(100.0%) | Tablets |       |        |         |