

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | pm98770       | Batch NO. :          | SHM2501210            |
| Print Interval : | 50000 Tablets | Running Time :       | 22 Hour 15 Min 25 Sec |
| Total Counter :  | 2601960       | Total Running Time : | 1034 Hour 21 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 26       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.500    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 21.71    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 20.87    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 17.03    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 16.70    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 16.36    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 12.52    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 11.69    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |               |       |            |       |       |            |    |
|----------------------|---------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.46          | kN    | Min :      | 4.32  | kN    | Punch NO.; | 70 |
| SD :                 | 1.20          | %     | Max :      | 4.62  | kN    | Punch NO.; | 75 |
| Mean Main Pressure : | 16.91         | kN    | Min :      | 15.43 | kN    | Punch NO.; | 70 |
| SD :                 | 3.30          | %     | Max :      | 18.48 | kN    | Punch NO.; | 61 |
| Filling Depth :      | Increase ; 71 | Times | Decrease ; | 135   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60   | mm  | Pre Thickness ;  | 3.10     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.53   | mm  | Main Thickness ; | 2.03     | mm           |
| Filling Depth :               | 8.19   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.757  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 22     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 25     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 13 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 36 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 14 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 2601960             | Tablets |       |        |         |
| A.W.C Counter : | 2596786             | Tablets |       |        |         |
| HEP :           | 26 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 2596760<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Sertra 100mgSTSH4B00  |
| User ID :        | pm98770       | Batch NO. :          | SHM2501210            |
| Print Interval : | 50000 Tablets | Running Time :       | 22 Hour 15 Min 25 Sec |
| Total Counter :  | 2601960       | Total Running Time : | 1034 Hour 21 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 25       | %       | Manual ;           | 48           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 3000000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.500    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 30       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 22.00    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 21.15    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 17.26    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 16.92    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 16.58    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 12.69    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 11.85    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 4.43           | kN    | Min :      | 4.35  | kN    | Punch NO.; | 14 |
| SD :                 | 1.20           | %     | Max :      | 4.61  | kN    | Punch NO.; | 36 |
| Mean Main Pressure : | 16.74          | kN    | Min :      | 15.47 | kN    | Punch NO.; | 17 |
| SD :                 | 3.70           | %     | Max :      | 18.67 | kN    | Punch NO.; | 75 |
| Filling Depth :      | Increase ; 226 | Times | Decrease ; | 234   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 7        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60   | mm  | Pre Thickness ;  | 3.10     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.64   | mm  | Main Thickness ; | 2.14     | mm           |
| Filling Depth :               | 8.02   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.757  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 22     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 25     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.20 | kN |             |    |
|                  | Max ;     | 0.31 | kN | Max Punch ; | 13 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 36 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 14 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 2601960             | Tablets |       |        |         |
| A.W.C Counter : | 2596786             | Tablets |       |        |         |
| HEP :           | 26 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 2596760<br>(100.0%) | Tablets |       |        |         |