

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pm98770       | Batch NO. :          | EWM2507750             |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 16 Min 25 Sec  |
| Total Counter :  | 2189522       | Total Running Time : | 1071 Hour 10 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 50       | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 4500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 30       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 10       | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 28.08    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 27.00    | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 22.03    | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 21.60    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 21.16    | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 16.20    | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 15.12    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 4.78       | kN | Min : | 4.48       | kN | Punch NO.; | 33 |
| SD :                 | 2.10       | %  | Max : | 5.00       | kN | Punch NO.; | 45 |
| Mean Main Pressure : | 21.64      | kN | Min : | 19.70      | kN | Punch NO.; | 46 |
| SD :                 | 3.90       | %  | Max : | 23.89      | kN | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 50 | Times | Decrease ; | 57 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 14       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.35    | mm  | Pre Thickness ;  | 3.36     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.40    | mm  | Main Thickness ; | 2.41     | mm           |
| Filling Depth :               | 7.82    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.829   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 11      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 36 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 2  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 18 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 2189522             | Tablets |       |        |         |
| A.W.C Counter : | 2181100             | Tablets |       |        |         |
| HEP :           | 29 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 2181071<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pm98770       | Batch NO. :          | EWM2507750             |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 16 Min 25 Sec  |
| Total Counter :  | 2189522       | Total Running Time : | 1071 Hour 10 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 50        | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm              |
|                                    |        | 8.5-14mm  |         |                    | Max : | 12.00 mm             |
| Target Quantity :                  |        | 4500000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.600     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 30        | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 10        | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 5.0       | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 27.77     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 26.70     | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 21.79     | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 21.36     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 20.93     | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 16.02     | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 14.95     | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |           |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 4.71       | kN | Min : | 4.40       | kN | Punch NO.; | 37 |
| SD :                 | 2.20       | %  | Max : | 5.02       | kN | Punch NO.; | 1  |
| Mean Main Pressure : | 21.32      | kN | Min : | 19.82      | kN | Punch NO.; | 46 |
| SD :                 | 4.00       | %  | Max : | 23.77      | kN | Punch NO.; | 13 |
| Filling Depth :      | Increase ; | 33 | Times | Decrease ; | 37 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 14       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.30    | mm  | Pre Thickness ;  | 3.31     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.42    | mm  | Main Thickness ; | 2.43     | mm           |
| Filling Depth :               | 7.55    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.829   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 11      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.19 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 36 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 2  |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 2189522             | Tablets |       |        |         |
| A.W.C Counter : | 2181100             | Tablets |       |        |         |
| HEP :           | 29 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 2181071<br>(100.0%) | Tablets |       |        |         |