

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2507950             |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 40 Min 05 Sec  |
| Total Counter :  | 1073218       | Total Running Time : | 1116 Hour 54 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 45                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 45                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 18                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 10.92             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 10.50             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 8.57              | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 8.40              | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 8.23              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 6.30              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 5.88              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |            |      |       |            |    |
|----------------------|--------------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | 3.12         | kN    | Min :      | 2.94 | kN    | Punch NO.; | 21 |
| SD :                 | 2.80         | %     | Max :      | 3.44 | kN    | Punch NO.; | 8  |
| Mean Main Pressure : | 8.35         | kN    | Min :      | 7.30 | kN    | Punch NO.; | 41 |
| SD :                 | 5.50         | %     | Max :      | 9.55 | kN    | Punch NO.; | 21 |
| Filling Depth :      | Increase ; 3 | Times | Decrease ; | 13   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.44    | mm  | Pre Thickness ;  | 3.44     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.25    | mm  | Main Thickness ; | 2.25     | mm           |
| Filling Depth :               | 6.87    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.700   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 31      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 6       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 15 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 6  |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1073218             | Tablets |       |           |         |
| A.W.C Counter : | 1060423             | Tablets |       |           |         |
| HEP :           | 75 (0.0%)           | Tablets | LEP : | 31 (0.0%) | Tablets |
| GOOD :          | 1060317<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2507950             |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 40 Min 05 Sec  |
| Total Counter :  | 1073218       | Total Running Time : | 1116 Hour 54 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 38                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 45                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 18                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 10.59             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.18             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 8.31              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.15              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 7.98              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.11              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 5.70              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.76         | kN    | Min :        | 2.57  | kN | Punch NO.; | 13 |
| SD :                 | 3.80         | %     | Max :        | 3.13  | kN | Punch NO.; | 8  |
| Mean Main Pressure : | 7.95         | kN    | Min :        | 7.11  | kN | Punch NO.; | 3  |
| SD :                 | 5.10         | %     | Max :        | 8.93  | kN | Punch NO.; | 21 |
| Filling Depth :      | Increase ; 6 | Times | Decrease ; 6 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.42    | mm  | Pre Thickness ;  | 3.42     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.46    | mm  | Main Thickness ; | 2.46     | mm           |
| Filling Depth :               | 6.47    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.700   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 31      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 6       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 15 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 52 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 1073218             | Tablets |       |           |         |
| A.W.C Counter : | 1060423             | Tablets |       |           |         |
| HEP :           | 75 (0.0%)           | Tablets | LEP : | 31 (0.0%) | Tablets |
| GOOD :          | 1060317<br>(100.0%) | Tablets |       |           |         |