

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2507950             |
| Print Interval : | 50000 Tablets | Running Time :       | 17 Hour 09 Min 05 Sec  |
| Total Counter :  | 2397780       | Total Running Time : | 1126 Hour 22 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 45                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 100               | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 11.43             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 10.99             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 8.97              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 8.79              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 8.61              | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 6.59              | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 6.15              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 3.29       | kN | Min : | 3.11       | kN | Punch NO.; | 26 |
| SD :                 | 3.30       | %  | Max : | 3.57       | kN | Punch NO.; | 11 |
| Mean Main Pressure : | 8.69       | kN | Min : | 7.48       | kN | Punch NO.; | 42 |
| SD :                 | 7.20       | %  | Max : | 10.64      | kN | Punch NO.; | 21 |
| Filling Depth :      | Increase ; | 34 | Times | Decrease ; | 38 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.44    | mm  | Pre Thickness ;  | 3.44     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.25    | mm  | Main Thickness ; | 2.25     | mm           |
| Filling Depth :               | 6.85    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.819   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 27      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 88      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 2397780             | Tablets |       |           |         |
| A.W.C Counter : | 2382801             | Tablets |       |           |         |
| HEP :           | 285 (0.0%)          | Tablets | LEP : | 31 (0.0%) | Tablets |
| GOOD :          | 2382485<br>(100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2507950             |
| Print Interval : | 50000 Tablets | Running Time :       | 17 Hour 09 Min 05 Sec  |
| Total Counter :  | 2397780       | Total Running Time : | 1126 Hour 22 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 38                | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 100               | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 10.90             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 10.48             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 8.55              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 8.39              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 8.22              | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 6.29              | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 5.87              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 3.15       | kN | Min : | 2.93       | kN | Punch NO.; | 30 |
| SD :                 | 3.60       | %  | Max : | 3.46       | kN | Punch NO.; | 53 |
| Mean Main Pressure : | 8.29       | kN | Min : | 7.05       | kN | Punch NO.; | 44 |
| SD :                 | 6.40       | %  | Max : | 9.46       | kN | Punch NO.; | 53 |
| Filling Depth :      | Increase ; | 22 | Times | Decrease ; | 23 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.35    | mm  | Pre Thickness ;  | 3.35     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.46    | mm  | Main Thickness ; | 2.46     | mm           |
| Filling Depth :               | 6.47    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.819   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 27      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 88      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.03 | kN | Max Punch ; | 15 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 53 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 2397780             | Tablets |       |           |         |
| A.W.C Counter : | 2382854             | Tablets |       |           |         |
| HEP :           | 285 (0.0%)          | Tablets | LEP : | 31 (0.0%) | Tablets |
| GOOD :          | 2382538<br>(100.0%) | Tablets |       |           |         |