

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2507950             |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 29 Min 18 Sec  |
| Total Counter :  | 3398294       | Total Running Time : | 1133 Hour 42 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 45                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 100               | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 10.39             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 9.99              | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 8.15              | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 7.99              | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 7.83              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 6.00              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 5.60              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 3.41          | kN    | Min :         | 3.16  | kN | Punch NO.; | 43 |
| SD :                 | 3.00          | %     | Max :         | 3.67  | kN | Punch NO.; | 9  |
| Mean Main Pressure : | 7.91          | kN    | Min :         | 6.58  | kN | Punch NO.; | 3  |
| SD :                 | 7.30          | %     | Max :         | 9.60  | kN | Punch NO.; | 14 |
| Filling Depth :      | Increase ; 92 | Times | Decrease ; 97 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.40    | mm  | Pre Thickness ;  | 3.40     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.30    | mm  | Main Thickness ; | 2.30     | mm           |
| Filling Depth :               | 6.85    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.845   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 46      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 67      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 52 |

### Tablet Data (N/A)

|                 |                  |         |       |           |         |
|-----------------|------------------|---------|-------|-----------|---------|
| Total Counter : | 3398294          | Tablets |       |           |         |
| A.W.C Counter : | 3383105          | Tablets |       |           |         |
| HEP :           | 493 (0.0%)       | Tablets | LEP : | 68 (0.0%) | Tablets |
| GOOD :          | 3382544 (100.0%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2507950             |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 29 Min 18 Sec  |
| Total Counter :  | 3398294       | Total Running Time : | 1133 Hour 42 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 38                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 100               | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 10.54             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.13             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 8.27              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.11              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 7.95              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.08              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 5.68              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 3.14       | kN | Min : | 2.86       | kN | Punch NO.; | 30 |
| SD :                 | 3.60       | %  | Max : | 3.44       | kN | Punch NO.; | 36 |
| Mean Main Pressure : | 8.14       | kN | Min : | 6.88       | kN | Punch NO.; | 3  |
| SD :                 | 6.70       | %  | Max : | 9.49       | kN | Punch NO.; | 14 |
| Filling Depth :      | Increase ; | 65 | Times | Decrease ; | 67 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 127200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.35    | mm  | Pre Thickness ;  | 3.35     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.46    | mm  | Main Thickness ; | 2.46     | mm           |
| Filling Depth :               | 6.43    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.845   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 46      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 67      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 26 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 52 |

### Tablet Data (N/A)

|                 |                     |         |       |           |         |
|-----------------|---------------------|---------|-------|-----------|---------|
| Total Counter : | 3398294             | Tablets |       |           |         |
| A.W.C Counter : | 3383105             | Tablets |       |           |         |
| HEP :           | 493 (0.0%)          | Tablets | LEP : | 68 (0.0%) | Tablets |
| GOOD :          | 3382544<br>(100.0%) | Tablets |       |           |         |