

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2508350             |
| Print Interval : | 50000 Tablets | Running Time :       | 15 Hour 42 Min 17 Sec  |
| Total Counter :  | 2801318       | Total Running Time : | 1183 Hour 15 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 69                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 100               | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 10.80             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 10.39             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 8.48              | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 8.31              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 8.14              | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 6.23              | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 5.82              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |      |       |            |      |       |            |    |
|----------------------|------------|------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.04 | kN    | Min :      | 2.90 | kN    | Punch NO.; | 13 |
|                      |            | 1.80 | %     | Max :      | 3.18 | kN    | Punch NO.; | 17 |
| Mean Main Pressure : |            | 8.30 | kN    | Min :      | 7.11 | kN    | Punch NO.; | 46 |
|                      | SD :       | 6.10 | %     | Max :      | 9.81 | kN    | Punch NO.; | 21 |
| Filling Depth :      | Increase ; | 1303 | Times | Decrease ; | 1370 | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 3.00   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.47   | mm  | Pre Thickness ;  | 3.47     | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.26   | mm  | Main Thickness ; | 2.26     | mm           |
| Filling Depth :               | 7.30   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.850  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 43     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 24     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 37 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                 |         |       |            |         |
|-----------------|-----------------|---------|-------|------------|---------|
| Total Counter : | 2801318         | Tablets |       |            |         |
| A.W.C Counter : | 2793872         | Tablets |       |            |         |
| HEP :           | 1566 (0.1%)     | Tablets | LEP : | 501 (0.0%) | Tablets |
| GOOD :          | 2791805 (99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2508350             |
| Print Interval : | 50000 Tablets | Running Time :       | 15 Hour 42 Min 17 Sec  |
| Total Counter :  | 2801318       | Total Running Time : | 1183 Hour 15 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 68                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 60                | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 100               | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 10.71             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.30             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %  | 8.40              | kN      | Yes 10             | Times        |           |
| Ref :                              |        | 8.24              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.07              | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 % | 6.18              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 % | 5.77              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.97       | kN  | Min : | 2.76       | kN  | Punch NO.; | 13 |
| SD :                 | 2.80       | %   | Max : | 3.18       | kN  | Punch NO.; | 52 |
| Mean Main Pressure : | 8.19       | kN  | Min : | 7.22       | kN  | Punch NO.; | 3  |
| SD :                 | 5.90       | %   | Max : | 9.39       | kN  | Punch NO.; | 21 |
| Filling Depth :      | Increase ; | 789 | Times | Decrease ; | 794 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.39   | mm  | Pre Thickness ;  | 3.40     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.46   | mm  | Main Thickness ; | 2.47     | mm           |
| Filling Depth :               | 6.94   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.850  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 43     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 24     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 37 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 13 |

### Tablet Data (N/A)

|                 |                 |         |       |            |         |
|-----------------|-----------------|---------|-------|------------|---------|
| Total Counter : | 2801318         | Tablets |       |            |         |
| A.W.C Counter : | 2793872         | Tablets |       |            |         |
| HEP :           | 1566 (0.1%)     | Tablets | LEP : | 501 (0.0%) | Tablets |
| GOOD :          | 2791805 (99.9%) | Tablets |       |            |         |