

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2508950             |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 20 Min 28 Sec  |
| Total Counter :  | 1485988       | Total Running Time : | 1201 Hour 3 Min        |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 61                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam C<br>8.5-14mm |         | Cam Range :        | Min : | 6.50 mm              |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max : | 12.00 mm             |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 120               | Min     | Supply Time ;      | 1.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 120               | Min     | Supply Time ;      | 1.8   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 27.99             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 26.91             | kN      | Yes                | 10    | Rotations 10 Tablets |
| HCP :                              | 2.0 %  | 21.96             | kN      | Yes                | 10    | Times                |
| Ref :                              |        | 21.53             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 21.10             | kN      | Yes                | 10    | Times                |
| LEP :                              | 25.0 % | 16.15             | kN      | Yes                | 10    | Rotations 10 Tablets |
| LSP :                              | 30.0 % | 15.07             | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.75       | kN | Min : | 2.58       | kN | Punch NO.; | 13 |
| SD :                 | 2.30       | %  | Max : | 2.88       | kN | Punch NO.; | 36 |
| Mean Main Pressure : | 21.27      | kN | Min : | 19.94      | kN | Punch NO.; | 3  |
| SD :                 | 3.00       | %  | Max : | 22.95      | kN | Punch NO.; | 21 |
| Filling Depth :      | Increase ; | 78 | Times | Decrease ; | 83 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 17       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.26    | mm  | Main Thickness ; | 2.26     | mm           |
| Filling Depth :               | 7.55    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.835   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 70      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 71      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 30 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 43 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1485988             | Tablets |       |        |         |
| A.W.C Counter : | 1479570             | Tablets |       |        |         |
| HEP :           | 10 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1479560<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2508950             |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 20 Min 28 Sec  |
| Total Counter :  | 1485988       | Total Running Time : | 1201 Hour 3 Min        |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 71                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 12.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 1.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 1.8          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 26.55             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 25.53             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 20.83             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 20.43             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 20.02             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 15.32             | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 14.30             | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.58       | kN | Min : | 2.31       | kN | Punch NO.; | 13 |
| SD :                 | 4.90       | %  | Max : | 2.93       | kN | Punch NO.; | 51 |
| Mean Main Pressure : | 19.27      | kN | Min : | 17.07      | kN | Punch NO.; | 45 |
| SD :                 | 5.30       | %  | Max : | 22.20      | kN | Punch NO.; | 21 |
| Filling Depth :      | Increase ; | 75 | Times | Decrease ; | 64 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.39    | mm  | Pre Thickness ;  | 3.40     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.46    | mm  | Main Thickness ; | 2.47     | mm           |
| Filling Depth :               | 7.39    | mm  | Current Cam :    | Cam C    |              |
|                               |         |     |                  | 8.5-14mm |              |
| Main Air pressure :           | 0.835   | MPa |                  |          |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 70      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 71      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 30 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 7  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 43 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1485988             | Tablets |       |        |         |
| A.W.C Counter : | 1479570             | Tablets |       |        |         |
| HEP :           | 10 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1479560<br>(100.0%) | Tablets |       |        |         |