

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 57 Min 59 Sec  |
| Total Counter :  | 350240        | Total Running Time : | 1220 Hour 53 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 78                | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 1.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 1.8          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 14.26             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 13.71             | kN      | Yes 10             | Rotations 10 | Tablets   |
| HCP :                              | 2.0 %       | 11.19             | kN      | Yes 10             | Times        |           |
| Ref :                              |             | 10.97             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.75             | kN      | Yes 10             | Times        |           |
| LEP :                              | 25.0 %      | 8.23              | kN      | Yes 10             | Rotations 10 | Tablets   |
| LSP :                              | 30.0 %      | 7.68              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.18       | kN  | Min : | 2.96       | kN  | Punch NO.; | 22 |
| SD :                 | 4.50       | %   | Max : | 3.87       | kN  | Punch NO.; | 50 |
| Mean Main Pressure : | 11.39      | kN  | Min : | 9.98       | kN  | Punch NO.; | 41 |
| SD :                 | 6.00       | %   | Max : | 13.67      | kN  | Punch NO.; | 10 |
| Filling Depth :      | Increase ; | 345 | Times | Decrease ; | 299 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 22                | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.34    | mm  | Main Thickness ; | 2.35              | mm           |
| Filling Depth :               | 7.34    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.812   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 350240            | Tablets |       |            |         |
| A.W.C Counter : | 346764            | Tablets |       |            |         |
| HEP :           | 162 (0.0%)        | Tablets | LEP : | 298 (0.1%) | Tablets |
| GOOD :          | 346304<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | sj97904       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 57 Min 59 Sec  |
| Total Counter :  | 350240        | Total Running Time : | 1220 Hour 53 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 78                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 1.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 1.8          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 13.18             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 12.67             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 10.34             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 10.14             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 9.94              | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 7.60              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 7.10              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.91       | kN  | Min : | 2.68       | kN  | Punch NO.; | 27 |
| SD :                 | 4.50       | %   | Max : | 3.29       | kN  | Punch NO.; | 42 |
| Mean Main Pressure : | 10.40      | kN  | Min : | 9.06       | kN  | Punch NO.; | 30 |
| SD :                 | 5.90       | %   | Max : | 11.83      | kN  | Punch NO.; | 20 |
| Filling Depth :      | Increase ; | 212 | Times | Decrease ; | 217 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 22       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.40    | mm  | Pre Thickness ;  | 3.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.48    | mm  | Main Thickness ; | 2.49     | mm           |
| Filling Depth :               | 6.90    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.812   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.27 | kN | Max Punch ; | 29 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.24 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 350240            | Tablets |       |            |         |
| A.W.C Counter : | 346764            | Tablets |       |            |         |
| HEP :           | 162 (0.0%)        | Tablets | LEP : | 298 (0.1%) | Tablets |
| GOOD :          | 346304<br>(99.9%) | Tablets |       |            |         |