

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | am98755       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 44 Min 29 Sec  |
| Total Counter :  | 666560        | Total Running Time : | 1222 Hour 39 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 71                | %       | Manual ;           | 15           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 1.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 1.8          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 13.44             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 12.92             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 10.54             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 10.34             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 10.13             | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 7.75              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 7.24              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | None        |                   |         |                    |              |            |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.00  | kN    | Min :      | 2.80  | kN    | Punch NO.; | 17 |
|                      |            | 3.30  | %     | Max :      | 3.27  | kN    | Punch NO.; | 20 |
| Mean Main Pressure : |            | 10.35 | kN    | Min :      | 8.52  | kN    | Punch NO.; | 29 |
|                      | SD :       | 6.20  | %     | Max :      | 12.02 | kN    | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 84    | Times | Decrease ; | 84    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.34    | mm  | Main Thickness ; | 2.35     | mm           |
| Filling Depth :               | 7.03    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.799   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.33 | kN |             |    |
|                  | Max ;     | 0.42 | kN | Max Punch ; | 24 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 666560            | Tablets |       |            |         |
| A.W.C Counter : | 661849            | Tablets |       |            |         |
| HEP :           | 363 (0.1%)        | Tablets | LEP : | 441 (0.1%) | Tablets |
| GOOD :          | 661045<br>(99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | am98755       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 44 Min 29 Sec  |
| Total Counter :  | 666560        | Total Running Time : | 1222 Hour 39 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 70                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 120               | Min     | Supply Time ;      | 1.5          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 120               | Min     | Supply Time ;      | 1.8          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 12.74             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 12.25             | kN      | Yes                | 10 Rotations | 10 Tablets |
| HCP :                              | 2.0 %       | 10.00             | kN      | Yes                | 10 Times     |            |
| Ref :                              |             | 9.80              | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 9.61              | kN      | Yes                | 10 Times     |            |
| LEP :                              | 25.0 %      | 7.35              | kN      | Yes                | 10 Rotations | 10 Tablets |
| LSP :                              | 30.0 %      | 6.86              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations  |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.78       | kN | Min : | 2.51       | kN | Punch NO.; | 17 |
| SD :                 | 4.20       | %  | Max : | 3.08       | kN | Punch NO.; | 46 |
| Mean Main Pressure : | 9.76       | kN | Min : | 8.36       | kN | Punch NO.; | 29 |
| SD :                 | 6.20       | %  | Max : | 10.98      | kN | Punch NO.; | 26 |
| Filling Depth :      | Increase ; | 36 | Times | Decrease ; | 26 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.40    | mm  | Pre Thickness ;  | 3.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.48    | mm  | Main Thickness ; | 2.49     | mm           |
| Filling Depth :               | 6.79    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.799   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.18 | kN |             |    |
|                  | Max ;     | 0.26 | kN | Max Punch ; | 44 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.33 | kN |             |    |
|                  | Max ;     | 0.42 | kN | Max Punch ; | 24 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 666560            | Tablets |       |            |         |
| A.W.C Counter : | 661902            | Tablets |       |            |         |
| HEP :           | 363 (0.1%)        | Tablets | LEP : | 441 (0.1%) | Tablets |
| GOOD :          | 661098<br>(99.9%) | Tablets |       |            |         |