

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | ms88737       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 20 Hour 34 Min 46 Sec  |
| Total Counter :  | 3669570       | Total Running Time : | 1239 Hour 30 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 72                | %       | Manual ;           | 15    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 4500000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 20                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 14.28             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 13.73             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 11.20             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 10.98             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 10.76             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 8.24              | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 7.69              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.09  | kN    | Min :      | 2.88  | kN    | Punch NO.; | 22 |
|                      |            | 2.90  | %     | Max :      | 3.38  | kN    | Punch NO.; | 38 |
| Mean Main Pressure : |            | 10.68 | kN    | Min :      | 9.16  | kN    | Punch NO.; | 15 |
|                      | SD :       | 6.40  | %     | Max :      | 12.20 | kN    | Punch NO.; | 20 |
| Filling Depth :      | Increase ; | 341   | Times | Decrease ; | 334   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.46    | mm  | Pre Thickness ;  | 3.46     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.34    | mm  | Main Thickness ; | 2.35     | mm           |
| Filling Depth :               | 7.18    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.806   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.16 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 10 |

### Tablet Data (N/A)

|                 |                 |         |       |            |         |
|-----------------|-----------------|---------|-------|------------|---------|
| Total Counter : | 3669570         | Tablets |       |            |         |
| A.W.C Counter : | 3655399         | Tablets |       |            |         |
| HEP :           | 1243 (0.0%)     | Tablets | LEP : | 999 (0.0%) | Tablets |
| GOOD :          | 3653157 (99.9%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | ms88737       | Batch NO. :          | EWM2509450             |
| Print Interval : | 50000 Tablets | Running Time :       | 20 Hour 34 Min 46 Sec  |
| Total Counter :  | 3669570       | Total Running Time : | 1239 Hour 30 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |            |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|------------|
| Feeder :                           | Auto ;      | 64                | %       | Manual ;           | 20           | RPM        |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm    |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 8.00 mm    |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |            |
| Oil Lubrication :                  |             |                   |         |                    |              |            |
| Upper Punch Interval Time(S1) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec        |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec        |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |            |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |            |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |            |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |            |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |            |
|                                    |             |                   |         | Stop Condition ;   |              |            |
| HSP :                              | 30.0 %      | 13.58             | kN      | No                 |              |            |
| HEP :                              | 25.0 %      | 13.06             | kN      | Yes                | 10 Rotations | 15 Tablets |
| HCP :                              | 2.0 %       | 10.65             | kN      | Yes                | 15 Times     |            |
| Ref :                              |             | 10.45             | kN      |                    |              |            |
| LCP :                              | 2.0 %       | 10.24             | kN      | Yes                | 15 Times     |            |
| LEP :                              | 25.0 %      | 7.83              | kN      | Yes                | 10 Rotations | 15 Tablets |
| LSP :                              | 30.0 %      | 7.31              | kN      | No                 |              |            |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |            |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |            |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2 Rotations  |            |
| Empty Punch No. :                  | N           |                   |         |                    |              |            |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 2.95       | kN  | Min : | 2.72       | kN  | Punch NO.; | 27 |
| SD :                 | 3.00       | %   | Max : | 3.16       | kN  | Punch NO.; | 41 |
| Mean Main Pressure : | 10.18      | kN  | Min : | 8.60       | kN  | Punch NO.; | 30 |
| SD :                 | 6.40       | %   | Max : | 11.57      | kN  | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 231 | Times | Decrease ; | 226 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.40    | mm  | Pre Thickness ;  | 3.41     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.50    | mm  | Main Thickness ; | 2.51     | mm           |
| Filling Depth :               | 6.91    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.806   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.16 | kN |             |    |
|                  | Max ;     | 0.25 | kN | Max Punch ; | 42 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 27 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 10 |

### Tablet Data (N/A)

|                 |                    |         |       |            |         |
|-----------------|--------------------|---------|-------|------------|---------|
| Total Counter : | 3669570            | Tablets |       |            |         |
| A.W.C Counter : | 3655399            | Tablets |       |            |         |
| HEP :           | 1243 (0.0%)        | Tablets | LEP : | 999 (0.0%) | Tablets |
| GOOD :          | 3653157<br>(99.9%) | Tablets |       |            |         |