

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2512250             |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 52 Min 18 Sec  |
| Total Counter :  | 1578692       | Total Running Time : | 1509 Hour 40 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 105               | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70                | Min     | Supply Time ;      | 0.1          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60                | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 27.62             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 26.56             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 21.67             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 21.25             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 20.82             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 15.94             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 14.87             | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.86         | kN    | Min :        | 2.80  | kN | Punch NO.; | 30 |
| SD :                 | 0.80         | %     | Max :        | 2.94  | kN | Punch NO.; | 11 |
| Mean Main Pressure : | 21.42        | kN    | Min :        | 20.13 | kN | Punch NO.; | 34 |
| SD :                 | 2.70         | %     | Max :        | 22.82 | kN | Punch NO.; | 11 |
| Filling Depth :      | Increase ; 0 | Times | Decrease ; 7 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 29                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.60    | mm  | Main Thickness ; | 2.61              | mm           |
| Filling Depth :               | 7.79    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.797   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 20      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 34      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 23 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 4  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1578692             | Tablets |       |            |         |
| A.W.C Counter : | 1571073             | Tablets |       |            |         |
| HEP :           | 61 (0.0%)           | Tablets | LEP : | 629 (0.0%) | Tablets |
| GOOD :          | 1570383<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | pb97604       | Batch NO. :          | EWM2512250             |
| Print Interval : | 50000 Tablets | Running Time :       | 08 Hour 52 Min 18 Sec  |
| Total Counter :  | 1578692       | Total Running Time : | 1509 Hour 40 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 100       | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 4500000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 70        | Min     | Supply Time ;      | 0.1          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 60        | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |        | 30.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 31.66     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 30.44     | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 24.84     | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 24.36     | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 23.87     | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 18.27     | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 17.05     | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.64       | kN | Min : | 2.58       | kN | Punch NO.; | 2  |
| SD :                 | 1.20       | %  | Max : | 2.71       | kN | Punch NO.; | 10 |
| Mean Main Pressure : | 24.77      | kN | Min : | 22.31      | kN | Punch NO.; | 34 |
| SD :                 | 3.60       | %  | Max : | 26.38      | kN | Punch NO.; | 13 |
| Filling Depth :      | Increase ; | 0  | Times | Decrease ; | 1  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 28                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.60    | mm  | Main Thickness ; | 2.61              | mm           |
| Filling Depth :               | 7.53    | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.797   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 20      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 34      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 16 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 4  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1578692             | Tablets |       |            |         |
| A.W.C Counter : | 1571073             | Tablets |       |            |         |
| HEP :           | 61 (0.0%)           | Tablets | LEP : | 629 (0.0%) | Tablets |
| GOOD :          | 1570383<br>(100.0%) | Tablets |       |            |         |