

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | yu98249       | Batch NO. :          | EWM2512250             |
| Print Interval : | 50000 Tablets | Running Time :       | 18 Hour 58 Min 41 Sec  |
| Total Counter :  | 3379784       | Total Running Time : | 1519 Hour 47 Min       |
| Station :        | Station 1     |                      |                        |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 105               | %       | Manual ;           | 15           | RPM       |
| Filling Cam :                      |             | Cam C<br>8.5-14mm |         | Cam Range :        | Min :        | 6.50 mm   |
| Target Quantity :                  |             | 4500000           | Tablets |                    | Max :        | 12.00 mm  |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70                | Min     | Supply Time ;      | 0.2          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 60                | Min     | Supply Time ;      | 0.2          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 28.69             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 27.59             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 22.51             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 22.07             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 21.63             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 16.55             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 15.45             | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.89          | kN    | Min :         | 2.83  | kN | Punch NO.; | 34 |
| SD :                 | 0.80          | %     | Max :         | 2.94  | kN | Punch NO.; | 11 |
| Mean Main Pressure : | 21.89         | kN    | Min :         | 20.53 | kN | Punch NO.; | 34 |
| SD :                 | 2.60          | %     | Max :         | 23.04 | kN | Punch NO.; | 26 |
| Filling Depth :      | Increase ; 37 | Times | Decrease ; 56 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 29       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.60    | mm  | Main Thickness ; | 2.61     | mm           |
| Filling Depth :               | 7.87    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.787   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 19      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 23 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 3379784             | Tablets |       |            |         |
| A.W.C Counter : | 3371906             | Tablets |       |            |         |
| HEP :           | 135 (0.0%)          | Tablets | LEP : | 630 (0.0%) | Tablets |
| GOOD :          | 3371141<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                        |
|------------------|---------------|----------------------|------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 500<br>MG,STEW1A00 |
| User ID :        | yu98249       | Batch NO. :          | EWM2512250             |
| Print Interval : | 50000 Tablets | Running Time :       | 18 Hour 58 Min 41 Sec  |
| Total Counter :  | 3379784       | Total Running Time : | 1519 Hour 47 Min       |
| Station :        | Station 2     |                      |                        |

### Setting Value

|                                    |        |           |         |                    |       |                      |
|------------------------------------|--------|-----------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 100       | %       | Manual ;           | 20    | RPM                  |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min : | 6.50 mm              |
|                                    |        | 8.5-14mm  |         |                    | Max : | 12.00 mm             |
| Target Quantity :                  |        | 4500000   | Tablets |                    |       |                      |
| Air Pressure Low Limit :           |        | 0.400     | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |           |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 70        | Min     | Supply Time ;      | 0.2   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 60        | Min     | Supply Time ;      | 0.2   | Sec                  |
| Powder Supply Time :               |        | 30.0      | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1         | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00      | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |       |                      |
|                                    |        |           |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 30.89     | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 29.70     | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 24.24     | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 23.76     | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 23.29     | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 17.82     | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 16.63     | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |           |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00     | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |           |         |                    |       |                      |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.66  | kN    | Min :      | 2.58  | kN    | Punch NO.; | 34 |
|                      |            | 1.10  | %     | Max :      | 2.73  | kN    | Punch NO.; | 27 |
| Mean Main Pressure : | SD :       | 23.65 | kN    | Min :      | 21.15 | kN    | Punch NO.; | 34 |
|                      |            | 3.30  | %     | Max :      | 25.06 | kN    | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 36    | Times | Decrease ; | 42    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28      | RPM | Capacity :       | 178080   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 28       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.65    | mm  | Main Thickness ; | 2.66     | mm           |
| Filling Depth :               | 7.46    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.787   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 19      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 13      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 17 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 23 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 3379784             | Tablets |       |            |         |
| A.W.C Counter : | 3371906             | Tablets |       |            |         |
| HEP :           | 135 (0.0%)          | Tablets | LEP : | 630 (0.0%) | Tablets |
| GOOD :          | 3371141<br>(100.0%) | Tablets |       |            |         |