

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	SERTRA50MG,STJL1A00
User ID :	ks99394	Batch NO. :	JLM2500450
Print Interval :	50000 Tablets	Running Time :	00 Hour 50 Min 03 Sec
Total Counter :	155884	Total Running Time :	1552 Hour 10 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	50	%	Manual ;	35	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		700000	Tablets			
Air Pressure Low Limit :		0.500	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		60	Min	Supply Time ;	0.2	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	9.08	kN	No		
HEP :	25.0 %	8.73	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	7.12	kN	Yes 15	Times	
Ref :		6.98	kN			
LCP :	2.0 %	6.84	kN	Yes 15	Times	
LEP :	25.0 %	5.24	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	4.89	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.62	kN	Min :	2.59	kN	Punch NO.;	61
SD :	0.50	%	Max :	2.64	kN	Punch NO.;	36
Mean Main Pressure :	7.00	kN	Min :	6.58	kN	Punch NO.;	7
SD :	2.90	%	Max :	7.50	kN	Punch NO.;	39
Filling Depth :	Increase ; 36	Times	Decrease ; 36	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	30	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.49	mm	Pre Thickness ;	3.50	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.70	mm	Main Thickness ;	1.71	mm
Filling Depth :	5.90	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.806	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	44	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.71	kN		
	Max ;	0.92	kN	Max Punch ;	65
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.10	kN	Max Punch ;	3
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.08	kN	Max Punch ;	74

Tablet Data (N/A)

Total Counter :	155884	Tablets			
A.W.C Counter :	134931	Tablets			
HEP :	915 (0.7%)	Tablets	LEP :	64 (0.0%)	Tablets
GOOD :	133952 (99.3%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	SERTRA50MG,STJL1A00
User ID :	ks99394	Batch NO. :	JLM2500450
Print Interval :	50000 Tablets	Running Time :	00 Hour 50 Min 03 Sec
Total Counter :	155884	Total Running Time :	1552 Hour 10 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	50	%	Manual ;	16	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		700000	Tablets			
Air Pressure Low Limit :		0.500	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	0.2	Sec
Lower Punch Interval Time(S2) ;		60	Min	Supply Time ;	0.2	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	7.23	kN	No		
HEP :	24.0 %	6.90	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	5.67	kN	Yes 15	Times	
Ref :		5.56	kN			
LCP :	2.0 %	5.45	kN	Yes 15	Times	
LEP :	24.0 %	4.23	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	3.89	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.39	kN	Min :	2.36	kN	Punch NO.;	17
SD :	0.70	%	Max :	2.47	kN	Punch NO.;	36
Mean Main Pressure :	5.45	kN	Min :	4.77	kN	Punch NO.;	14
SD :	5.90	%	Max :	6.69	kN	Punch NO.;	36
Filling Depth :	Increase ;	79	Times	Decrease ;	67	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	30	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	15	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	6.49	mm	Pre Thickness ;	3.50	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	4.75	mm	Main Thickness ;	1.76	mm
Filling Depth :	6.11	mm	Current Cam :	Cam B 4.5-10mm	
Main Air pressure :	0.806	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	44	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.73	kN		
	Max ;	0.92	kN	Max Punch ;	65
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.10	kN	Max Punch ;	3
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.08	kN	Max Punch ;	74

Tablet Data (N/A)

Total Counter :	155884	Tablets			
A.W.C Counter :	134931	Tablets			
HEP :	915 (0.7%)	Tablets	LEP :	64 (0.0%)	Tablets
GOOD :	133952 (99.3%)	Tablets			