

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA100MG<br>STJL1B00 |
| User ID :        | rp98927       | Batch NO. :          | JLM2500410              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 38 Min 55 Sec   |
| Total Counter :  | 83274         | Total Running Time : | 1560 Hour 53 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 44                | %       | Manual ;           | 33           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 650000            | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 45                | Min     | Supply Time ;      | 0.6          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 14.56             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 14.00             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 11.42             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 11.20             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 10.98             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 8.40              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 7.84              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |               |       |    |            |    |
|----------------------|--------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.92         | kN    | Min :         | 2.87  | kN | Punch NO.; | 3  |
| SD :                 | 0.70         | %     | Max :         | 2.98  | kN | Punch NO.; | 76 |
| Mean Main Pressure : | 11.01        | kN    | Min :         | 10.01 | kN | Punch NO.; | 15 |
| SD :                 | 3.40         | %     | Max :         | 11.82 | kN | Punch NO.; | 40 |
| Filling Depth :      | Increase ; 7 | Times | Decrease ; 11 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 30      | RPM | Capacity :       | 277200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 13       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.97    | mm  | Main Thickness ; | 1.97     | mm           |
| Filling Depth :               | 7.47    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.798   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 15      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.49 | kN |             |    |
|                  | Max ;     | 0.63 | kN | Max Punch ; | 33 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 42 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 83274             | Tablets |       |        |         |
| A.W.C Counter : | 67488             | Tablets |       |        |         |
| HEP :           | 2 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 67486<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA100MG<br>STJL1B00 |
| User ID :        | rp98927       | Batch NO. :          | JLM2500410              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 38 Min 55 Sec   |
| Total Counter :  | 83274         | Total Running Time : | 1560 Hour 53 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 52                | %       | Manual ;           | 16           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 650000            | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 45                | Min     | Supply Time ;      | 0.6          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 11.40             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 10.97             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 8.95              | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 8.77              | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 8.60              | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 6.58              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 6.14              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |              |       |    |            |    |
|----------------------|--------------|-------|--------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.70         | kN    | Min :        | 2.65  | kN | Punch NO.; | 13 |
| SD :                 | 0.90         | %     | Max :        | 2.77  | kN | Punch NO.; | 7  |
| Mean Main Pressure : | 8.71         | kN    | Min :        | 7.98  | kN | Punch NO.; | 16 |
| SD :                 | 4.70         | %     | Max :        | 9.75  | kN | Punch NO.; | 40 |
| Filling Depth :      | Increase ; 1 | Times | Decrease ; 7 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 30      | RPM | Capacity :       | 277200   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 16       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.18    | mm  | Main Thickness ; | 2.19     | mm           |
| Filling Depth :               | 7.54    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.798   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 15      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.49 | kN |             |    |
|                  | Max ;     | 0.63 | kN | Max Punch ; | 33 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 42 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 83274             | Tablets |       |        |         |
| A.W.C Counter : | 67488             | Tablets |       |        |         |
| HEP :           | 2 (0.0%)          | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 67486<br>(100.0%) | Tablets |       |        |         |