

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA100MG<br>STJL1B00 |
| User ID :        | sj97904       | Batch NO. :          | JLM2500410              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 44 Min 58 Sec   |
| Total Counter :  | 375590        | Total Running Time : | 1561 Hour 59 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 48                | %       | Manual ;           | 33           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 650000            | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 3.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 15.44             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 14.84             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 12.11             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 11.88             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 11.64             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 8.91              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 8.31              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.94       | kN | Min : | 2.91       | kN | Punch NO.; | 39 |
| SD :                 | 0.40       | %  | Max : | 2.98       | kN | Punch NO.; | 36 |
| Mean Main Pressure : | 11.97      | kN | Min : | 11.34      | kN | Punch NO.; | 22 |
| SD :                 | 2.50       | %  | Max : | 13.29      | kN | Punch NO.; | 36 |
| Filling Depth :      | Increase ; | 44 | Times | Decrease ; | 65 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.97    | mm  | Main Thickness ; | 1.97     | mm           |
| Filling Depth :               | 7.25    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.806   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 3       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.42 | kN |             |    |
|                  | Max ;     | 0.73 | kN | Max Punch ; | 38 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 47 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |           |         |
|-----------------|-------------------|---------|-------|-----------|---------|
| Total Counter : | 375590            | Tablets |       |           |         |
| A.W.C Counter : | 358671            | Tablets |       |           |         |
| HEP :           | 212 (0.1%)        | Tablets | LEP : | 49 (0.0%) | Tablets |
| GOOD :          | 358410<br>(99.9%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA100MG<br>STJL1B00 |
| User ID :        | sj97904       | Batch NO. :          | JLM2500410              |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 44 Min 58 Sec   |
| Total Counter :  | 375590        | Total Running Time : | 1561 Hour 59 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 50                | %       | Manual ;           | 16    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 650000            | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 25                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 3.0   | Sec                  |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 12.17             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 11.71             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 9.55              | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 9.36              | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 9.18              | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 7.02              | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 6.56              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.71       | kN | Min : | 2.68       | kN | Punch NO.; | 21 |
| SD :                 | 0.80       | %  | Max : | 2.81       | kN | Punch NO.; | 75 |
| Mean Main Pressure : | 9.49       | kN | Min : | 8.97       | kN | Punch NO.; | 14 |
| SD :                 | 3.80       | %  | Max : | 10.91      | kN | Punch NO.; | 36 |
| Filling Depth :      | Increase ; | 18 | Times | Decrease ; | 21 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 231000            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.00    | mm  | Pre Thickness ;  | 4.01              | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.18    | mm  | Main Thickness ; | 2.19              | mm           |
| Filling Depth :               | 7.70    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.806   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 3       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | On                |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.42 | kN |             |    |
|                  | Max ;     | 0.72 | kN | Max Punch ; | 38 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 47 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                   |         |       |           |         |
|-----------------|-------------------|---------|-------|-----------|---------|
| Total Counter : | 375590            | Tablets |       |           |         |
| A.W.C Counter : | 358671            | Tablets |       |           |         |
| HEP :           | 212 (0.1%)        | Tablets | LEP : | 49 (0.0%) | Tablets |
| GOOD :          | 358410<br>(99.9%) | Tablets |       |           |         |