

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 20 MG,STKQ1A00
User ID :	sj97904	Batch NO. :	KQM2501920
Print Interval :	50000 Tablets	Running Time :	01 Hour 19 Min 09 Sec
Total Counter :	300170	Total Running Time :	1579 Hour 11 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	19	%	Manual ;	18	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		5500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		25	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	14.71	kN	No		
HEP :	25.0 %	14.15	kN	Yes	10	Rotations 15 Tablets
HCP :	2.0 %	11.55	kN	Yes	15	Times
Ref :		11.32	kN			
LCP :	2.0 %	11.09	kN	Yes	15	Times
LEP :	25.0 %	8.49	kN	Yes	10	Rotations 15 Tablets
LSP :	30.0 %	7.92	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.01	kN	Min :	1.56	kN	Punch NO.;	49
SD :	8.50	%	Max :	2.42	kN	Punch NO.;	40
Mean Main Pressure :	11.35	kN	Min :	8.82	kN	Punch NO.;	33
SD :	7.80	%	Max :	13.18	kN	Punch NO.;	2
Filling Depth :	Increase ;	391	Times	Decrease ;	378	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	5	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.10	mm	Pre Thickness ;	1.10	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.55	mm	Main Thickness ;	0.55	mm
Filling Depth :	6.67	mm	Current Cam :	Cam B	
Main Air pressure :	0.737	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.09	kN		
	Max ;	0.14	kN	Max Punch ;	28
Lower Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Ejection Force :	Average ;	0.17	kN		
	Max ;	0.33	kN	Max Punch ;	60

Tablet Data (N/A)

Total Counter :	300170	Tablets			
A.W.C Counter :	272779	Tablets			
HEP :	1783 (0.7%)	Tablets	LEP :	1475 (0.5%)	Tablets
GOOD :	269521 (98.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	PANTO 20 MG,STKQ1A00
User ID :	sj97904	Batch NO. :	KQM2501920
Print Interval :	50000 Tablets	Running Time :	01 Hour 19 Min 09 Sec
Total Counter :	300170	Total Running Time :	1579 Hour 11 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	13	%	Manual ;	20	RPM
Filling Cam :		Cam B 4.5-10mm		Cam Range :	Min :	2.50 mm
Target Quantity :		5500000	Tablets		Max :	8.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		25	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	16.12	kN	No		
HEP :	25.0 %	15.50	kN	Yes	10 Rotations	15 Tablets
HCP :	2.0 %	12.64	kN	Yes	15 Times	
Ref :		12.40	kN			
LCP :	2.0 %	12.15	kN	Yes	15 Times	
LEP :	25.0 %	9.30	kN	Yes	10 Rotations	15 Tablets
LSP :	30.0 %	8.68	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.76	kN	Min :	1.39	kN	Punch NO.;	11
SD :	9.20	%	Max :	2.21	kN	Punch NO.;	40
Mean Main Pressure :	12.68	kN	Min :	9.48	kN	Punch NO.;	33
SD :	8.50	%	Max :	14.81	kN	Punch NO.;	39
Filling Depth :	Increase ;	220	Times	Decrease ;	222	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	3	RPM
Pre Upper Punch Position ;	3.00	mm			
Pre Lower Punch Position ;	4.10	mm	Pre Thickness ;	1.10	mm
Main Upper Punch Position ;	3.00	mm			
Main Lower Punch Position ;	3.55	mm	Main Thickness ;	0.55	mm
Filling Depth :	6.54	mm	Current Cam :	Cam B	
Main Air pressure :	0.737	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.09	kN		
	Max ;	0.14	kN	Max Punch ;	28
Lower Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Ejection Force :	Average ;	0.17	kN		
	Max ;	0.33	kN	Max Punch ;	60

Tablet Data (N/A)

Total Counter :	300170	Tablets			
A.W.C Counter :	272779	Tablets			
HEP :	1783 (0.7%)	Tablets	LEP :	1475 (0.5%)	Tablets
GOOD :	269521 (98.8%)	Tablets			