

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2501920              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 33 Min 50 Sec   |
| Total Counter :  | 700366        | Total Running Time : | 1581 Hour 26 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 15                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 25                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 15.71             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 15.10             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 12.32             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 12.08             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 11.84             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 9.06              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 8.46              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.11           | kN    | Min :      | 1.74  | kN    | Punch NO.; | 21 |
| SD :                 | 8.70           | %     | Max :      | 2.56  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : | 12.13          | kN    | Min :      | 10.21 | kN    | Punch NO.; | 60 |
| SD :                 | 9.00           | %     | Max :      | 15.16 | kN    | Punch NO.; | 34 |
| Filling Depth :      | Increase ; 807 | Times | Decrease ; | 769   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 30      | RPM | Capacity :       | 234000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.10    | mm  | Pre Thickness ;  | 1.10     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.55    | mm  | Main Thickness ; | 0.55     | mm           |
| Filling Depth :               | 6.27    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.743   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 6       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Ejection Force : | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 45 |

### Tablet Data (N/A)

|                 |                   |         |       |             |         |
|-----------------|-------------------|---------|-------|-------------|---------|
| Total Counter : | 700366            | Tablets |       |             |         |
| A.W.C Counter : | 661260            | Tablets |       |             |         |
| HEP :           | 3600 (0.5%)       | Tablets | LEP : | 4814 (0.7%) | Tablets |
| GOOD :          | 652846<br>(98.7%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | pb97604       | Batch NO. :          | KQM2501920              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 33 Min 50 Sec   |
| Total Counter :  | 700366        | Total Running Time : | 1581 Hour 26 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 15                | %       | Manual ;           | 14           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 25                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 15.62             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 15.02             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 12.26             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 12.02             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 11.78             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 9.01              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 8.41              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.65  | kN    | Min :      | 1.31  | kN    | Punch NO.; | 13 |
|                      |            | 11.30 | %     | Max :      | 2.03  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : |            | 12.29 | kN    | Min :      | 9.33  | kN    | Punch NO.; | 7  |
|                      | SD :       | 8.80  | %     | Max :      | 14.92 | kN    | Punch NO.; | 34 |
| Filling Depth :      | Increase ; | 617   | Times | Decrease ; | 603   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 30      | RPM | Capacity :       | 234000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.10    | mm  | Pre Thickness ;  | 1.10     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.55    | mm  | Main Thickness ; | 0.55     | mm           |
| Filling Depth :               | 6.42    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.743   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 6       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 7       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 41 |
| Lower Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Ejection Force : | Average ; | 0.14 | kN |             |    |
|                  | Max ;     | 0.22 | kN | Max Punch ; | 45 |

### Tablet Data (N/A)

|                 |                   |         |       |             |         |
|-----------------|-------------------|---------|-------|-------------|---------|
| Total Counter : | 700366            | Tablets |       |             |         |
| A.W.C Counter : | 661260            | Tablets |       |             |         |
| HEP :           | 3600 (0.5%)       | Tablets | LEP : | 4814 (0.7%) | Tablets |
| GOOD :          | 652846<br>(98.7%) | Tablets |       |             |         |