

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	SERTRA100MG STSH4B00
User ID :	sj97904	Batch NO. :	SHM2502310
Print Interval :	50000 Tablets	Running Time :	00 Hour 11 Min 37 Sec
Total Counter :	50014	Total Running Time :	1647 Hour 58 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	36	%	Manual ;	20	RPM
Filling Cam :		Cam C 8.5-14mm		Cam Range :	Min :	6.50 mm
Target Quantity :		3000000	Tablets		Max :	12.00 mm
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	18.51	kN	No		
HEP :	30.0 %	18.51	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	14.52	kN	Yes 15	Times	
Ref :		14.24	kN			
LCP :	2.0 %	13.95	kN	Yes 15	Times	
LEP :	30.0 %	9.97	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	9.97	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	1	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.44	kN	Min :	1.40	kN	Punch NO.;	14
SD :	1.00	%	Max :	1.47	kN	Punch NO.;	1
Mean Main Pressure :	14.29	kN	Min :	13.25	kN	Punch NO.;	71
SD :	3.10	%	Max :	15.08	kN	Punch NO.;	34
Filling Depth :	Increase ;	428	Times	Decrease ;	522	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	28	RPM	Capacity :	258720	Tablets/hour
Feeder : Status ;	Auto		Speed ;	10	RPM
Pre Upper Punch Position ;	2.50	mm			
Pre Lower Punch Position ;	5.58	mm	Pre Thickness ;	3.08	mm
Main Upper Punch Position ;	2.50	mm			
Main Lower Punch Position ;	4.50	mm	Main Thickness ;	2.00	mm
Filling Depth :	9.82	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.801	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	23	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	On	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.01	kN		
	Max ;	0.05	kN	Max Punch ;	37
Lower Punch :	Average ;	0.10	kN		
	Max ;	0.11	kN	Max Punch ;	30
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	49

Tablet Data (N/A)

Total Counter :	50014	Tablets			
A.W.C Counter :	44250	Tablets			
HEP :	0 (0%)	Tablets	LEP :	69 (0.2%)	Tablets
GOOD :	44181 (99.8%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	SERTRA100MG STSH4B00
User ID :	sj97904	Batch NO. :	SHM2502310
Print Interval :	50000 Tablets	Running Time :	00 Hour 11 Min 37 Sec
Total Counter :	50014	Total Running Time :	1647 Hour 58 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	28	%	Manual ;	14	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		3000000	Tablets			
Air Pressure Low Limit :		0.600	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		60	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		30.0	Sec			
Initial Reject Time :		1	Tablets			
Upper Punch Tightness High Limit :		2.00	kN			
Lower Punch Tightness High Limit :		4.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	17.06	kN	No		
HEP :	30.0 %	17.06	kN	Yes 10	Rotations 15	Tablets
HCP :	2.0 %	13.39	kN	Yes 15	Times	
Ref :		13.12	kN			
LCP :	2.0 %	12.86	kN	Yes 15	Times	
LEP :	30.0 %	9.19	kN	Yes 10	Rotations 15	Tablets
LSP :	30.0 %	9.19	kN	No		
SD limit :	35.0 %			Yes		
Pre HSP		30.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	1	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	1.00	kN	Min :	0.95	kN	Punch NO.;	42
SD :	2.10	%	Max :	1.06	kN	Punch NO.;	4
Mean Main Pressure :	13.07	kN	Min :	12.09	kN	Punch NO.;	29
SD :	3.80	%	Max :	14.09	kN	Punch NO.;	8
Filling Depth :	Increase ;	466	Times	Decrease ;	427	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	28	RPM	Capacity :	258720	Tablets/hour
Feeder : Status ;	Auto		Speed ;	8	RPM
Pre Upper Punch Position ;	2.50	mm			
Pre Lower Punch Position ;	5.50	mm	Pre Thickness ;	3.00	mm
Main Upper Punch Position ;	2.50	mm			
Main Lower Punch Position ;	4.53	mm	Main Thickness ;	2.03	mm
Filling Depth :	9.60	mm	Current Cam :	Cam C	
				8.5-14mm	
Main Air pressure :	0.801	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	23	Min			
Lower Punch Remain Time (S2):	5	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	On	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.01	kN		
	Max ;	0.05	kN	Max Punch ;	37
Lower Punch :	Average ;	0.10	kN		
	Max ;	0.11	kN	Max Punch ;	72
Ejection Force :	Average ;	0.05	kN		
	Max ;	0.08	kN	Max Punch ;	49

Tablet Data (N/A)

Total Counter :	50014	Tablets			
A.W.C Counter :	44250	Tablets			
HEP :	0 (0%)	Tablets	LEP :	69 (0.2%)	Tablets
GOOD :	44181 (99.8%)	Tablets			