

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | sj97904       | Batch NO. :          | KQM2502120              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 08 Min 40 Sec   |
| Total Counter :  | 500248        | Total Running Time : | 1696 Hour 7 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |           |
|------------------------------------|--------|-------------------|---------|--------------------|-------|-----------|
| Feeder :                           | Auto ; | 25                | %       | Manual ;           | 20    | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm   |
| Target Quantity :                  |        | 5500000           | Tablets |                    | Max : | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |       |           |
| Oil Lubrication :                  |        |                   |         |                    |       |           |
| Upper Punch Interval Time(S1) ;    |        | 58                | Min     | Supply Time ;      | 1.0   | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 1.6   | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |       |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |       |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |       |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |           |
|                                    |        |                   |         | Stop Condition ;   |       |           |
| HSP :                              | 30.0 % | 17.00             | kN      | No                 |       |           |
| HEP :                              | 25.0 % | 16.35             | kN      | Yes                | 10    | Rotations |
| HCP :                              | 2.0 %  | 13.34             | kN      | Yes                | 15    | Times     |
| Ref :                              |        | 13.08             | kN      |                    |       |           |
| LCP :                              | 2.0 %  | 12.82             | kN      | Yes                | 15    | Times     |
| LEP :                              | 25.0 % | 9.81              | kN      | Yes                | 10    | Rotations |
| LSP :                              | 30.0 % | 9.16              | kN      | No                 |       |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations |
| Empty Punch No. :                  | None   |                   |         |                    |       |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.69  | kN    | Min :      | 2.21  | kN    | Punch NO.; | 38 |
|                      |            | 7.40  | %     | Max :      | 3.11  | kN    | Punch NO.; | 25 |
| Mean Main Pressure : |            | 13.85 | kN    | Min :      | 10.71 | kN    | Punch NO.; | 17 |
|                      | SD :       | 7.40  | %     | Max :      | 15.77 | kN    | Punch NO.; | 25 |
| Filling Depth :      | Increase ; | 246   | Times | Decrease ; | 244   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 30      | RPM | Capacity :       | 234000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.35    | mm  | Pre Thickness ;  | 1.35     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 0.60     | mm           |
| Filling Depth :               | 6.42    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.804   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 49      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.21 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 37 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 9  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 500248            | Tablets |       |            |         |
| A.W.C Counter : | 491687            | Tablets |       |            |         |
| HEP :           | 895 (0.2%)        | Tablets | LEP : | 452 (0.1%) | Tablets |
| GOOD :          | 490340<br>(99.7%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | sj97904       | Batch NO. :          | KQM2502120              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 08 Min 40 Sec   |
| Total Counter :  | 500248        | Total Running Time : | 1696 Hour 7 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 29                | %       | Manual ;           | 14           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 58                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 1.6          | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 17.48             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 16.81             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 13.72             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 13.45             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 13.18             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 10.09             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 9.41              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.96       | kN  | Min : | 1.51       | kN  | Punch NO.; | 17 |
| SD :                 | 8.80       | %   | Max : | 2.43       | kN  | Punch NO.; | 42 |
| Mean Main Pressure : | 13.28      | kN  | Min : | 9.92       | kN  | Punch NO.; | 17 |
| SD :                 | 8.00       | %   | Max : | 16.13      | kN  | Punch NO.; | 42 |
| Filling Depth :      | Increase ; | 250 | Times | Decrease ; | 261 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 30      | RPM | Capacity :       | 234000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.36    | mm  | Pre Thickness ;  | 1.36     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.66    | mm  | Main Thickness ; | 0.66     | mm           |
| Filling Depth :               | 6.37    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.804   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 49      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 61 |
| Lower Punch :    | Average ; | 0.21 | kN |             |    |
|                  | Max ;     | 0.37 | kN | Max Punch ; | 37 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 9  |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 500248            | Tablets |       |            |         |
| A.W.C Counter : | 491687            | Tablets |       |            |         |
| HEP :           | 895 (0.2%)        | Tablets | LEP : | 452 (0.1%) | Tablets |
| GOOD :          | 490340<br>(99.7%) | Tablets |       |            |         |