

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | sj97904       | Batch NO. :          | KQM2502120              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 26 Min 07 Sec   |
| Total Counter :  | 5445366       | Total Running Time : | 1721 Hour 28 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 17                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 40                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 17.01             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.35             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 13.34             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 13.08             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 12.82             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 9.81              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 9.16              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.47       | kN  | Min : | 1.17       | kN  | Punch NO.; | 38 |
| SD :                 | 9.60       | %   | Max : | 1.88       | kN  | Punch NO.; | 57 |
| Mean Main Pressure : | 11.74      | kN  | Min : | 10.08      | kN  | Punch NO.; | 26 |
| SD :                 | 6.60       | %   | Max : | 14.19      | kN  | Punch NO.; | 57 |
| Filling Depth :      | Increase ; | 368 | Times | Decrease ; | 359 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 4                 | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.36   | mm  | Pre Thickness ;  | 1.37              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.60   | mm  | Main Thickness ; | 0.61              | mm           |
| Filling Depth :               | 6.71   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.719  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 4      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 2      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | On                |              |
| Initial Reject :              | On     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.06 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 16 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 54 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 54 |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 5445366            | Tablets |       |             |         |
| A.W.C Counter : | 5423559            | Tablets |       |             |         |
| HEP :           | 1877 (0.0%)        | Tablets | LEP : | 1909 (0.0%) | Tablets |
| GOOD :          | 5419773<br>(99.9%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 20<br>MG,STKQ1A00 |
| User ID :        | sj97904       | Batch NO. :          | KQM2502120              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 26 Min 07 Sec   |
| Total Counter :  | 5445366       | Total Running Time : | 1721 Hour 28 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 20                | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 5500000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.600             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 40                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 15                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 1                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 2.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 17.04             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 16.38             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 13.37             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 13.11             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 12.84             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 9.83              | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 9.17              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.07       | kN  | Min : | 0.90       | kN  | Punch NO.; | 38 |
| SD :                 | 13.70      | %   | Max : | 1.56       | kN  | Punch NO.; | 62 |
| Mean Main Pressure : | 8.24       | kN  | Min : | 6.43       | kN  | Punch NO.; | 17 |
| SD :                 | 15.10      | %   | Max : | 12.50      | kN  | Punch NO.; | 65 |
| Filling Depth :      | Increase ; | 564 | Times | Decrease ; | 488 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 5                 | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 4.36   | mm  | Pre Thickness ;  | 1.37              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.66   | mm  | Main Thickness ; | 0.67              | mm           |
| Filling Depth :               | 6.63   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.719  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 4      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 2      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | On                |              |
| Initial Reject :              | On     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 61 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 54 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.06 | kN | Max Punch ; | 54 |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 5445366            | Tablets |       |             |         |
| A.W.C Counter : | 5423559            | Tablets |       |             |         |
| HEP :           | 1877 (0.0%)        | Tablets | LEP : | 1909 (0.0%) | Tablets |
| GOOD :          | 5419773<br>(99.9%) | Tablets |       |             |         |