

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amio 200mg STAI3A00   |
| User ID :        | vg96371       | Batch NO. :          | AIM2500420            |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 10 Min 55 Sec |
| Total Counter :  | 89096         | Total Running Time : | 1750 Hour 27 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 122      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 26.14    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 25.13    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 20.51    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 20.11    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 19.70    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 15.08    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 14.07    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 3.19          | kN    | Min :         | 2.19  | kN | Punch NO.; | 13 |
| SD :                 | 11.20         | %     | Max :         | 4.01  | kN | Punch NO.; | 40 |
| Mean Main Pressure : | 19.90         | kN    | Min :         | 16.97 | kN | Punch NO.; | 13 |
| SD :                 | 5.30          | %     | Max :         | 21.99 | kN | Punch NO.; | 40 |
| Filling Depth :      | Increase ; 47 | Times | Decrease ; 49 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400             | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 18                | RPM          |
| Pre Upper Punch Position ;    | 2.20    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 3.40    | mm  | Pre Thickness ;  | 1.20              | mm           |
| Main Upper Punch Position ;   | 2.80    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.40    | mm  | Main Thickness ; | 0.60              | mm           |
| Filling Depth :               | 6.90    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.802   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 11      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.31 | kN |             |    |
|                  | Max ;     | 0.45 | kN | Max Punch ; | 51 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 43 |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 89096            | Tablets |       |            |         |
| A.W.C Counter : | 84083            | Tablets |       |            |         |
| HEP :           | 253 (0.3%)       | Tablets | LEP : | 119 (0.1%) | Tablets |
| GOOD :          | 83711<br>(99.6%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Amio 200mg STAI3A00   |
| User ID :        | vg96371       | Batch NO. :          | AIM2500420            |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 10 Min 55 Sec |
| Total Counter :  | 89096         | Total Running Time : | 1750 Hour 27 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 124      | %       | Manual ;           | 20           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 1500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.600    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 60       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 26.47    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 25.45    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 20.77    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 20.36    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 19.96    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 15.27    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 14.25    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.43          | kN    | Min :         | 1.78  | kN | Punch NO.; | 13 |
| SD :                 | 10.60         | %     | Max :         | 3.10  | kN | Punch NO.; | 53 |
| Mean Main Pressure : | 19.24         | kN    | Min :         | 17.08 | kN | Punch NO.; | 13 |
| SD :                 | 4.70          | %     | Max :         | 21.73 | kN | Punch NO.; | 53 |
| Filling Depth :      | Increase ; 43 | Times | Decrease ; 50 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 19       | RPM          |
| Pre Upper Punch Position ;    | 2.20    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 3.40    | mm  | Pre Thickness ;  | 1.20     | mm           |
| Main Upper Punch Position ;   | 2.80    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.40    | mm  | Main Thickness ; | 0.60     | mm           |
| Filling Depth :               | 6.88    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.802   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 11      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 2       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.31 | kN |             |    |
|                  | Max ;     | 0.47 | kN | Max Punch ; | 49 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 4  |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 43 |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 89096            | Tablets |       |            |         |
| A.W.C Counter : | 84083            | Tablets |       |            |         |
| HEP :           | 253 (0.3%)       | Tablets | LEP : | 119 (0.1%) | Tablets |
| GOOD :          | 83711<br>(99.6%) | Tablets |       |            |         |