

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40 MG<br>STKQ1B00 |
| User ID :        | os190259      | Batch NO. :          | KQM2507440              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 41 Min 14 Sec   |
| Total Counter :  | 111750        | Total Running Time : | 1798 Hour 51 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 55                | %       | Manual ;           | 10           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 35                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 34.82             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 33.48             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 27.32             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 26.79             | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 26.25             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 20.09             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 18.75             | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.22         | kN    | Min :      | 1.16  | kN    | Punch NO.; | 29 |
| SD :                 | 2.40         | %     | Max :      | 1.29  | kN    | Punch NO.; | 50 |
| Mean Main Pressure : | 26.78        | kN    | Min :      | 24.51 | kN    | Punch NO.; | 63 |
| SD :                 | 3.90         | %     | Max :      | 29.61 | kN    | Punch NO.; | 50 |
| Filling Depth :      | Increase ; 1 | Times | Decrease ; | 15    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600            | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.01              | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 3.95    | mm  | Main Thickness ; | 0.95              | mm           |
| Filling Depth :               | 4.82    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.796   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 10      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 13 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 64 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 34 |

### Tablet Data (N/A)

|                 |                  |         |       |           |         |
|-----------------|------------------|---------|-------|-----------|---------|
| Total Counter : | 111750           | Tablets |       |           |         |
| A.W.C Counter : | 94139            | Tablets |       |           |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 90 (0.1%) | Tablets |
| GOOD :          | 94049<br>(99.9%) | Tablets |       |           |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40 MG<br>STKQ1B00 |
| User ID :        | os190259      | Batch NO. :          | KQM2507440              |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 41 Min 14 Sec   |
| Total Counter :  | 111750        | Total Running Time : | 1798 Hour 51 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 55                | %       | Manual ;           | 10    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 35                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 33.13             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 31.86             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 25.99             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 25.48             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 24.97             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 19.11             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 17.84             | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | N      |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.16       | kN | Min : | 1.11       | kN | Punch NO.; | 18 |
| SD :                 | 1.90       | %  | Max : | 1.21       | kN | Punch NO.; | 3  |
| Mean Main Pressure : | 25.32      | kN | Min : | 23.83      | kN | Punch NO.; | 41 |
| SD :                 | 2.90       | %  | Max : | 27.37      | kN | Punch NO.; | 43 |
| Filling Depth :      | Increase ; | 1  | Times | Decrease ; | 1  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.00    | mm  | Pre Thickness ;  | 2.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.10    | mm  | Main Thickness ; | 1.10     | mm           |
| Filling Depth :               | 4.31    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.796   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 24      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 10      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 42 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 64 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 34 |

### Tablet Data (N/A)

|                 |                  |         |       |           |         |
|-----------------|------------------|---------|-------|-----------|---------|
| Total Counter : | 111750           | Tablets |       |           |         |
| A.W.C Counter : | 94139            | Tablets |       |           |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 90 (0.1%) | Tablets |
| GOOD :          | 94049<br>(99.9%) | Tablets |       |           |         |