

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40 MG<br>STKQ1B00 |
| User ID :        | cb83530       | Batch NO. :          | KQM2507540              |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 48 Min 30 Sec   |
| Total Counter :  | 1851044       | Total Running Time : | 1825 Hour 0 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 44                | %       | Manual ;           | 10    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 % | 34.06             | kN      | No                 |       |                      |
| HEP :                              | 25.0 % | 32.75             | kN      | Yes                | 10    | Rotations 15 Tablets |
| HCP :                              | 2.0 %  | 26.72             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 26.20             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 25.67             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 % | 19.65             | kN      | Yes                | 10    | Rotations 15 Tablets |
| LSP :                              | 30.0 % | 18.34             | kN      | No                 |       |                      |
| SD limit :                         | 35.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 1.47       | kN | Min : | 1.36       | kN | Punch NO.; | 16 |
| SD :                 | 3.50       | %  | Max : | 1.57       | kN | Punch NO.; | 49 |
| Mean Main Pressure : | 25.89      | kN | Min : | 23.63      | kN | Punch NO.; | 16 |
| SD :                 | 3.50       | %  | Max : | 27.95      | kN | Punch NO.; | 3  |
| Filling Depth :      | Increase ; | 25 | Times | Decrease ; | 33 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.30    | mm  | Pre Thickness ;  | 1.80     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.55    | mm  | Main Thickness ; | 1.05     | mm           |
| Filling Depth :               | 5.33    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.789   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 25      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.21 | kN |             |    |
|                  | Max ;     | 0.35 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 64 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 50 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1851044             | Tablets |       |            |         |
| A.W.C Counter : | 1840794             | Tablets |       |            |         |
| HEP :           | 53 (0.0%)           | Tablets | LEP : | 571 (0.0%) | Tablets |
| GOOD :          | 1840170<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO 40 MG<br>STKQ1B00 |
| User ID :        | cb83530       | Batch NO. :          | KQM2507540              |
| Print Interval : | 50000 Tablets | Running Time :       | 10 Hour 48 Min 30 Sec   |
| Total Counter :  | 1851044       | Total Running Time : | 1825 Hour 0 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 46                | %       | Manual ;           | 10           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2750000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 5.0               | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 4.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 36.02             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 34.63             | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %  | 28.26             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 27.71             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 27.15             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 20.78             | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 % | 19.40             | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.38  | kN    | Min :      | 1.30  | kN    | Punch NO.; | 40 |
|                      |            | 3.60  | %     | Max :      | 1.51  | kN    | Punch NO.; | 59 |
| Mean Main Pressure : | SD :       | 27.67 | kN    | Min :      | 25.39 | kN    | Punch NO.; | 40 |
|                      |            | 3.50  | %     | Max :      | 29.87 | kN    | Punch NO.; | 35 |
| Filling Depth :      | Increase ; | 53    | Times | Decrease ; | 52    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 10       | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.30    | mm  | Pre Thickness ;  | 1.80     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.60    | mm  | Main Thickness ; | 1.10     | mm           |
| Filling Depth :               | 4.43    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.789   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 25      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 3       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.22 | kN |             |    |
|                  | Max ;     | 0.35 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 64 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 50 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1851044             | Tablets |       |            |         |
| A.W.C Counter : | 1840794             | Tablets |       |            |         |
| HEP :           | 53 (0.0%)           | Tablets | LEP : | 571 (0.0%) | Tablets |
| GOOD :          | 1840170<br>(100.0%) | Tablets |       |            |         |