

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO40mgSTKQ1B00     |
| User ID :        | mr99799       | Batch NO. :          | KQM2507640            |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 26 Min 01 Sec |
| Total Counter :  | 1618374       | Total Running Time : | 1839 Hour 37 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 27       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 35.06    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 33.71    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 27.51    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 26.97    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 26.43    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 20.23    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 18.88    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.35  | kN    | Min :      | 1.21  | kN    | Punch NO.; | 55 |
|                      |            | 5.20  | %     | Max :      | 1.53  | kN    | Punch NO.; | 42 |
| Mean Main Pressure : |            | 27.33 | kN    | Min :      | 24.37 | kN    | Punch NO.; | 55 |
|                      | SD :       | 5.20  | %     | Max :      | 30.94 | kN    | Punch NO.; | 3  |
| Filling Depth :      | Increase ; | 312   | Times | Decrease ; | 323   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.30   | mm  | Pre Thickness ;  | 1.80     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.52   | mm  | Main Thickness ; | 1.02     | mm           |
| Filling Depth :               | 4.92   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.702  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 43     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 60 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 56 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 63 |

### Tablet Data (N/A)

|                 |                  |         |       |            |         |
|-----------------|------------------|---------|-------|------------|---------|
| Total Counter : | 1618374          | Tablets |       |            |         |
| A.W.C Counter : | 1613233          | Tablets |       |            |         |
| HEP :           | 103 (0.0%)       | Tablets | LEP : | 537 (0.0%) | Tablets |
| GOOD :          | 1612593 (100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO40mgSTKQ1B00     |
| User ID :        | mr99799       | Batch NO. :          | KQM2507640            |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 26 Min 01 Sec |
| Total Counter :  | 1618374       | Total Running Time : | 1839 Hour 37 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 38       | %       | Manual ;           | 10           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 38.19    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 36.72    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 29.97    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 29.38    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 28.79    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 22.03    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 20.56    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.32       | kN  | Min : | 1.20       | kN  | Punch NO.; | 12 |
| SD :                 | 4.80       | %   | Max : | 1.48       | kN  | Punch NO.; | 19 |
| Mean Main Pressure : | 29.37      | kN  | Min : | 27.32      | kN  | Punch NO.; | 12 |
| SD :                 | 4.10       | %   | Max : | 32.18      | kN  | Punch NO.; | 50 |
| Filling Depth :      | Increase ; | 212 | Times | Decrease ; | 227 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.30    | mm  | Pre Thickness ;  | 1.80     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.59    | mm  | Main Thickness ; | 1.09     | mm           |
| Filling Depth :               | 4.35    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.702   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 43      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 60 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 56 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 63 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1618374             | Tablets |       |            |         |
| A.W.C Counter : | 1613233             | Tablets |       |            |         |
| HEP :           | 103 (0.0%)          | Tablets | LEP : | 537 (0.0%) | Tablets |
| GOOD :          | 1612593<br>(100.0%) | Tablets |       |            |         |