

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO40mgSTKQ1B00     |
| User ID :        | pb97604       | Batch NO. :          | KQM2507740            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 17 Min 13 Sec |
| Total Counter :  | 1250700       | Total Running Time : | 1853 Hour 27 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 40       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 31.84    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 30.62    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 24.98    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 24.49    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 24.00    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 18.37    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 17.15    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.41           | kN    | Min :      | 1.24  | kN    | Punch NO.; | 24 |
| SD :                 | 3.80           | %     | Max :      | 1.53  | kN    | Punch NO.; | 39 |
| Mean Main Pressure : | 24.83          | kN    | Min :      | 22.60 | kN    | Punch NO.; | 24 |
| SD :                 | 3.90           | %     | Max :      | 27.30 | kN    | Punch NO.; | 43 |
| Filling Depth :      | Increase ; 672 | Times | Decrease ; | 685   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 9        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.30   | mm  | Pre Thickness ;  | 1.80     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.61   | mm  | Main Thickness ; | 1.11     | mm           |
| Filling Depth :               | 4.76   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.802  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 70     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 11 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 44 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 58 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1250700             | Tablets |       |        |         |
| A.W.C Counter : | 1243581             | Tablets |       |        |         |
| HEP :           | 56 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1243525<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO40mgSTKQ1B00     |
| User ID :        | pb97604       | Batch NO. :          | KQM2507740            |
| Print Interval : | 50000 Tablets | Running Time :       | 07 Hour 17 Min 13 Sec |
| Total Counter :  | 1250700       | Total Running Time : | 1853 Hour 27 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 28       | %       | Manual ;           | 10           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 5.0      | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 4.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 36.46    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 35.06    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 28.61    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 28.05    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 27.49    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 21.03    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 19.63    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.28           | kN    | Min :      | 1.18  | kN    | Punch NO.; | 56 |
| SD :                 | 4.30           | %     | Max :      | 1.47  | kN    | Punch NO.; | 43 |
| Mean Main Pressure : | 28.10          | kN    | Min :      | 26.01 | kN    | Punch NO.; | 56 |
| SD :                 | 3.50           | %     | Max :      | 31.40 | kN    | Punch NO.; | 43 |
| Filling Depth :      | Increase ; 360 | Times | Decrease ; | 364   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22     | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 2.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.30   | mm  | Pre Thickness ;  | 1.80     | mm           |
| Main Upper Punch Position ;   | 2.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.61   | mm  | Main Thickness ; | 1.11     | mm           |
| Filling Depth :               | 4.47   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.802  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 70     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 8      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 11 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 44 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 58 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1250700             | Tablets |       |        |         |
| A.W.C Counter : | 1243581             | Tablets |       |        |         |
| HEP :           | 56 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1243525<br>(100.0%) | Tablets |       |        |         |