

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO40mgSTKQ1B00     |
| User ID :        | rp98927       | Batch NO. :          | KQM2507840            |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 27 Min 49 Sec |
| Total Counter :  | 2138942       | Total Running Time : | 1874 Hour 39 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 36       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 32.56    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 31.31    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 25.55    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 25.05    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 24.55    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 18.79    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 17.53    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.05       | kN  | Min : | 1.00       | kN  | Punch NO.; | 58 |
| SD :                 | 2.00       | %   | Max : | 1.10       | kN  | Punch NO.; | 8  |
| Mean Main Pressure : | 25.28      | kN  | Min : | 23.21      | kN  | Punch NO.; | 58 |
| SD :                 | 3.80       | %   | Max : | 27.85      | kN  | Punch NO.; | 8  |
| Filling Depth :      | Increase ; | 111 | Times | Decrease ; | 129 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 8        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.85    | mm  | Pre Thickness ;  | 1.86     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.53    | mm  | Main Thickness ; | 1.03     | mm           |
| Filling Depth :               | 4.60    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.804   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 54      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 11 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 2138942             | Tablets |       |            |         |
| A.W.C Counter : | 2125950             | Tablets |       |            |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 574 (0.0%) | Tablets |
| GOOD :          | 2125374<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | PANTO40mgSTKQ1B00     |
| User ID :        | rp98927       | Batch NO. :          | KQM2507840            |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 27 Min 49 Sec |
| Total Counter :  | 2138942       | Total Running Time : | 1874 Hour 39 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 29       | %       | Manual ;           | 10           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 2750000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 32.89    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 31.62    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 25.80    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 25.30    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 24.79    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 18.97    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 17.71    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 1.04          | kN    | Min :         | 0.98  | kN | Punch NO.; | 17 |
| SD :                 | 2.60          | %     | Max :         | 1.11  | kN | Punch NO.; | 48 |
| Mean Main Pressure : | 25.29         | kN    | Min :         | 23.67 | kN | Punch NO.; | 31 |
| SD :                 | 3.50          | %     | Max :         | 27.13 | kN | Punch NO.; | 48 |
| Filling Depth :      | Increase ; 29 | Times | Decrease ; 31 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 171600   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 6        | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 4.85    | mm  | Pre Thickness ;  | 1.86     | mm           |
| Main Upper Punch Position ;   | 2.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 3.53    | mm  | Main Thickness ; | 1.03     | mm           |
| Filling Depth :               | 4.38    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.804   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 54      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 12      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |
| Lower Punch :    | Average ; | 0.13 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 26 |
| Ejection Force : | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.04 | kN | Max Punch ; | 63 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 2138942             | Tablets |       |            |         |
| A.W.C Counter : | 2125950             | Tablets |       |            |         |
| HEP :           | 2 (0.0%)            | Tablets | LEP : | 574 (0.0%) | Tablets |
| GOOD :          | 2125374<br>(100.0%) | Tablets |       |            |         |