

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET500MGSTMR7A00      |
| User ID :        | sj97904       | Batch NO. :          | MRM2502150            |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 55 Min 40 Sec |
| Total Counter :  | 1100530       | Total Running Time : | 1885 Hour 8 Min       |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 89       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 4500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 19.87    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 19.10    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 15.59    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 15.28    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.98    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 11.46    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.70    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.37  | kN    | Min :      | 1.30  | kN    | Punch NO.; | 8  |
|                      |            | 2.20  | %     | Max :      | 1.44  | kN    | Punch NO.; | 16 |
| Mean Main Pressure : |            | 15.34 | kN    | Min :      | 14.69 | kN    | Punch NO.; | 19 |
|                      | SD :       | 2.00  | %     | Max :      | 15.91 | kN    | Punch NO.; | 32 |
| Filling Depth :      | Increase ; | 54    | Times | Decrease ; | 70    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 22       | RPM          |
| Pre Upper Punch Position ;    | 2.49   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00   | mm  | Pre Thickness ;  | 3.51     | mm           |
| Main Upper Punch Position ;   | 2.49   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.60   | mm  | Main Thickness ; | 3.11     | mm           |
| Filling Depth :               | 7.82   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.792  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 47     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 10 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.30 | kN | Max Punch ; | 20 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1100530             | Tablets |       |            |         |
| A.W.C Counter : | 1088193             | Tablets |       |            |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 120 (0.0%) | Tablets |
| GOOD :          | 1088073<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET500MGSTMR7A00      |
| User ID :        | sj97904       | Batch NO. :          | MRM2502150            |
| Print Interval : | 50000 Tablets | Running Time :       | 06 Hour 55 Min 40 Sec |
| Total Counter :  | 1100530       | Total Running Time : | 1885 Hour 8 Min       |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 89       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 4500000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 20       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 30.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 19.38    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 18.63    | kN      | Yes 10             | Rotations 15 | Tablets   |
| HCP :                              | 2.0 %       | 15.21    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 14.91    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 14.61    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 11.18    | kN      | Yes 10             | Rotations 15 | Tablets   |
| LSP :                              | 30.0 %      | 10.44    | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.31  | kN    | Min :      | 1.21  | kN    | Punch NO.; | 39 |
|                      |            | 3.10  | %     | Max :      | 1.44  | kN    | Punch NO.; | 30 |
| Mean Main Pressure : |            | 14.87 | kN    | Min :      | 13.63 | kN    | Punch NO.; | 21 |
|                      | SD :       | 3.50  | %     | Max :      | 15.85 | kN    | Punch NO.; | 37 |
| Filling Depth :      | Increase ; | 39    | Times | Decrease ; | 46    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25     | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 22       | RPM          |
| Pre Upper Punch Position ;    | 2.49   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00   | mm  | Pre Thickness ;  | 3.51     | mm           |
| Main Upper Punch Position ;   | 2.49   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.56   | mm  | Main Thickness ; | 3.07     | mm           |
| Filling Depth :               | 7.40   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.792  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 47     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 5      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 10 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.30 | kN | Max Punch ; | 20 |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 1100530             | Tablets |       |            |         |
| A.W.C Counter : | 1088193             | Tablets |       |            |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 120 (0.0%) | Tablets |
| GOOD :          | 1088073<br>(100.0%) | Tablets |       |            |         |