

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA 50MG<br>STSH4A00 |
| User ID :        | sj97904       | Batch NO. :          | SHM2502350              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 57 Min 10 Sec   |
| Total Counter :  | 2501804       | Total Running Time : | 1944 Hour 40 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |       |                      |
|------------------------------------|-------------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ;      | 42                | %       | Manual ;           | 30    | RPM                  |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |       |                      |
| Oil Lubrication :                  |             |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.5   | Sec                  |
| Lower Punch Interval Time(S2) ;    |             | 20                | Min     | Supply Time ;      | 1.0   | Sec                  |
| Powder Supply Time :               |             | 30.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |             | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |       |                      |
|                                    |             |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 30.0 %      | 16.20             | kN      | No                 |       |                      |
| HEP :                              | 25.0 %      | 15.58             | kN      | Yes                | 10    | Rotations 20 Tablets |
| HCP :                              | 1.5 %       | 12.65             | kN      | Yes                | 15    | Times                |
| Ref :                              |             | 12.46             | kN      |                    |       |                      |
| LCP :                              | 1.5 %       | 12.28             | kN      | Yes                | 15    | Times                |
| LEP :                              | 25.0 %      | 9.35              | kN      | No                 | 35    | Rotations 35 Tablets |
| LSP :                              | 30.0 %      | 8.72              | kN      | No                 |       |                      |
| SD limit :                         | 35.0 %      |                   |         | Yes                |       |                      |
| Pre HSP                            |             | 30.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1     | Rotations            |
| Empty Punch No. :                  | None        |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |     |            |    |
|----------------------|------------|----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.97       | kN | Min : | 1.66       | kN  | Punch NO.; | 8  |
| SD :                 | 2.60       | %  | Max : | 2.05       | kN  | Punch NO.; | 70 |
| Mean Main Pressure : | 12.66      | kN | Min : | 9.24       | kN  | Punch NO.; | 8  |
| SD :                 | 4.80       | %  | Max : | 13.83      | kN  | Punch NO.; | 5  |
| Filling Depth :      | Increase ; | 64 | Times | Decrease ; | 101 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 12       | RPM          |
| Pre Upper Punch Position ;    | 3.00   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.60   | mm  | Pre Thickness ;  | 2.60     | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.11   | mm  | Main Thickness ; | 2.11     | mm           |
| Filling Depth :               | 5.78   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.769  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 21     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 63 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 8  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 2501804            | Tablets |       |             |         |
| A.W.C Counter : | 2482557            | Tablets |       |             |         |
| HEP :           | 178 (0.0%)         | Tablets | LEP : | 7157 (0.3%) | Tablets |
| GOOD :          | 2475222<br>(99.7%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA 50MG<br>STSH4A00 |
| User ID :        | sj97904       | Batch NO. :          | SHM2502350              |
| Print Interval : | 50000 Tablets | Running Time :       | 09 Hour 57 Min 10 Sec   |
| Total Counter :  | 2501804       | Total Running Time : | 1944 Hour 40 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 38                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 20                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 30.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 15.15             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 14.56             | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %  | 11.83             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 11.65             | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 11.48             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 8.74              | kN      | Yes 35             | Rotations 35 | Tablets   |
| LSP :                              | 30.0 % | 8.16              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |     |            |    |
|----------------------|------------|----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.93       | kN | Min : | 1.75       | kN  | Punch NO.; | 60 |
| SD :                 | 2.60       | %  | Max : | 2.02       | kN  | Punch NO.; | 38 |
| Mean Main Pressure : | 11.83      | kN | Min : | 10.03      | kN  | Punch NO.; | 60 |
| SD :                 | 4.20       | %  | Max : | 12.56      | kN  | Punch NO.; | 52 |
| Filling Depth :      | Increase ; | 93 | Times | Decrease ; | 141 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 258720   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 11       | RPM          |
| Pre Upper Punch Position ;    | 3.00   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.50     | mm           |
| Main Upper Punch Position ;   | 3.00   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.22   | mm  | Main Thickness ; | 2.22     | mm           |
| Filling Depth :               | 5.54   | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.769  | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 21     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 19     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | On       |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.18 | kN | Max Punch ; | 77 |
| Lower Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 8  |
| Ejection Force : | Average ; | 0.00 | kN |             |    |
|                  | Max ;     | 0.00 | kN | Max Punch ; | 0  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 2501804            | Tablets |       |             |         |
| A.W.C Counter : | 2482634            | Tablets |       |             |         |
| HEP :           | 178 (0.0%)         | Tablets | LEP : | 7157 (0.3%) | Tablets |
| GOOD :          | 2475299<br>(99.7%) | Tablets |       |             |         |