

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA 50MG<br>STSH4A00 |
| User ID :        | pr190041      | Batch NO. :          | SHM2502850              |
| Print Interval : | 50000 Tablets | Running Time :       | 23 Hour 02 Min 26 Sec   |
| Total Counter :  | 5653538       | Total Running Time : | 2073 Hour 0 Min         |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |                   |         |                    |              |           |
|------------------------------------|-------------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 50                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |             | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |             | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |             |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00              | kN      |                    |              |           |
|                                    |             |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 17.36             | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 16.69             | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %       | 13.55             | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 13.35             | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 13.15             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 10.01             | kN      | No 35              | Rotations 35 | Tablets   |
| LSP :                              | 30.0 %      | 9.35              | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |                   |         | Yes                |              |           |
| Pre HSP                            |             | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |                   |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.47       | kN  | Min : | 1.43       | kN  | Punch NO.; | 42 |
| SD :                 | 1.40       | %   | Max : | 1.53       | kN  | Punch NO.; | 23 |
| Mean Main Pressure : | 14.17      | kN  | Min : | 12.28      | kN  | Punch NO.; | 46 |
| SD :                 | 5.00       | %   | Max : | 15.88      | kN  | Punch NO.; | 62 |
| Filling Depth :      | Increase ; | 430 | Times | Decrease ; | 504 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 14                | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.27   | mm  | Main Thickness ; | 2.28              | mm           |
| Filling Depth :               | 5.91   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.748  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 40     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 21     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 51 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 18 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 5653538             | Tablets |       |            |         |
| A.W.C Counter : | 5647938             | Tablets |       |            |         |
| HEP :           | 5 (0.0%)            | Tablets | LEP : | 219 (0.0%) | Tablets |
| GOOD :          | 5647714<br>(100.0%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | SERTRA 50MG<br>STSH4A00 |
| User ID :        | pr190041      | Batch NO. :          | SHM2502850              |
| Print Interval : | 50000 Tablets | Running Time :       | 23 Hour 02 Min 26 Sec   |
| Total Counter :  | 5653538       | Total Running Time : | 2073 Hour 0 Min         |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 42                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 6000000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.400             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 17.49             | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 16.82             | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %  | 13.65             | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 13.45             | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 13.25             | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 10.09             | kN      | Yes 35             | Rotations 35 | Tablets   |
| LSP :                              | 30.0 % | 9.42              | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 30.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 1.17       | kN  | Min : | 1.14       | kN  | Punch NO.; | 44 |
| SD :                 | 1.40       | %   | Max : | 1.22       | kN  | Punch NO.; | 43 |
| Mean Main Pressure : | 13.92      | kN  | Min : | 10.62      | kN  | Punch NO.; | 63 |
| SD :                 | 6.10       | %   | Max : | 15.69      | kN  | Punch NO.; | 5  |
| Filling Depth :      | Increase ; | 786 | Times | Decrease ; | 804 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 28     | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 12                | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 5.50   | mm  | Pre Thickness ;  | 2.51              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 5.30   | mm  | Main Thickness ; | 2.31              | mm           |
| Filling Depth :               | 5.36   | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.748  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 40     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 21     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 54 |
| Lower Punch :    | Average ; | 0.08 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 18 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 33 |

### Tablet Data (N/A)

|                 |                     |         |       |            |         |
|-----------------|---------------------|---------|-------|------------|---------|
| Total Counter : | 5653538             | Tablets |       |            |         |
| A.W.C Counter : | 5647938             | Tablets |       |            |         |
| HEP :           | 5 (0.0%)            | Tablets | LEP : | 219 (0.0%) | Tablets |
| GOOD :          | 5647714<br>(100.0%) | Tablets |       |            |         |