

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 500MG<br>STHA1A00 |
| User ID :        | rp98927       | Batch NO. :          | HAM2501050              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 41 Min 02 Sec   |
| Total Counter :  | 350096        | Total Running Time : | 2077 Hour 40 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 98       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 750000   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 80       | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 13.39    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 12.88    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %       | 10.46    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 10.30    | kN      |                    |              |           |
| LCP :                              | 1.5 %       | 10.15    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 7.73     | kN      | No 35              | Rotations 35 | Tablets   |
| LSP :                              | 30.0 %      | 7.21     | kN      | No                 |              |           |
| SD limit :                         | 35.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 30.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 1 Rotations |          |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 0.97       | kN  | Min : | 0.95       | kN  | Punch NO.; | 12 |
| SD :                 | 1.20       | %   | Max : | 1.02       | kN  | Punch NO.; | 47 |
| Mean Main Pressure : | 10.08      | kN  | Min : | 7.61       | kN  | Punch NO.; | 21 |
| SD :                 | 10.50      | %   | Max : | 11.97      | kN  | Punch NO.; | 41 |
| Filling Depth :      | Increase ; | 624 | Times | Decrease ; | 561 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 15       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.50    | mm  | Pre Thickness ;  | 4.50     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.50    | mm  | Main Thickness ; | 2.50     | mm           |
| Filling Depth :               | 9.05    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.802   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 19      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 47 |
| Ejection Force : | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 22 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 350096            | Tablets |       |            |         |
| A.W.C Counter : | 324853            | Tablets |       |            |         |
| HEP :           | 540 (0.2%)        | Tablets | LEP : | 888 (0.3%) | Tablets |
| GOOD :          | 323425<br>(99.6%) | Tablets |       |            |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | METHO 500MG<br>STHA1A00 |
| User ID :        | rp98927       | Batch NO. :          | HAM2501050              |
| Print Interval : | 50000 Tablets | Running Time :       | 03 Hour 41 Min 02 Sec   |
| Total Counter :  | 350096        | Total Running Time : | 2077 Hour 40 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |           |         |                    |              |           |
|------------------------------------|--------|-----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 94        | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam C     |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |        | 8.5-14mm  |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |        | 750000    | Tablets |                    |              |           |
| Air Pressure Low Limit :           |        | 0.000     | MPa     |                    |              |           |
| Oil Lubrication :                  |        |           |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80        | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40        | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0      | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2         | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00      | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00      | kN      |                    |              |           |
|                                    |        |           |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 12.44     | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 11.96     | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 1.5 %  | 9.71      | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 9.57      | kN      |                    |              |           |
| LCP :                              | 1.5 %  | 9.43      | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 7.18      | kN      | Yes 35             | Rotations 35 | Tablets   |
| LSP :                              | 30.0 % | 6.70      | kN      | No                 |              |           |
| SD limit :                         | 35.0 % |           |         | Yes                |              |           |
| Pre HSP                            |        | 30.00     | kN      | Yes                |              |           |
| Mean Calculation range :           | 1      | Rotations |         | Calculation Pass : | 1            | Rotations |
| Empty Punch No. :                  | N      |           |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 0.74       | kN  | Min : | 0.72       | kN  | Punch NO.; | 25 |
| SD :                 | 2.00       | %   | Max : | 0.78       | kN  | Punch NO.; | 42 |
| Mean Main Pressure : | 9.20       | kN  | Min : | 7.99       | kN  | Punch NO.; | 30 |
| SD :                 | 6.70       | %   | Max : | 10.47      | kN  | Punch NO.; | 15 |
| Filling Depth :      | Increase ; | 511 | Times | Decrease ; | 507 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 15      | RPM | Capacity :       | 95400    | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 14       | RPM          |
| Pre Upper Punch Position ;    | 2.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.50    | mm  | Pre Thickness ;  | 4.50     | mm           |
| Main Upper Punch Position ;   | 2.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 4.66    | mm  | Main Thickness ; | 2.66     | mm           |
| Filling Depth :               | 8.11    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.802   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 19      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 1       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.07 | kN |             |    |
|                  | Max ;     | 0.12 | kN | Max Punch ; | 53 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 47 |
| Ejection Force : | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 41 |

### Tablet Data (N/A)

|                 |                   |         |       |            |         |
|-----------------|-------------------|---------|-------|------------|---------|
| Total Counter : | 350096            | Tablets |       |            |         |
| A.W.C Counter : | 324853            | Tablets |       |            |         |
| HEP :           | 540 (0.2%)        | Tablets | LEP : | 888 (0.3%) | Tablets |
| GOOD :          | 323425<br>(99.6%) | Tablets |       |            |         |