

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METHO 500MG STHA1A00
User ID :	kv190014	Batch NO. :	HAM2501150
Print Interval :	50000 Tablets	Running Time :	00 Hour 16 Min 28 Sec
Total Counter :	25792	Total Running Time :	2081 Hour 55 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	95	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		750000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	13.05	kN	No		
HEP :	25.0 %	12.55	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	10.24	kN	Yes 15	Times	
Ref :		10.04	kN			
LCP :	2.0 %	9.84	kN	Yes 15	Times	
LEP :	25.0 %	7.53	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	7.03	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	1.01	kN	Min :	0.97	kN	Punch NO.;	20
SD :	1.20	%	Max :	1.04	kN	Punch NO.;	29
Mean Main Pressure :	9.81	kN	Min :	8.59	kN	Punch NO.;	25
SD :	5.30	%	Max :	11.49	kN	Punch NO.;	15
Filling Depth :	Increase ;	1121	Times	Decrease ;	965	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.45	mm	Main Thickness ;	2.45	mm
Filling Depth :	7.63	mm	Current Cam :	Cam C	
Main Air pressure :	0.792	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	9	Min			
Lower Punch Remain Time (S2):	31	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.07	kN		
	Max ;	0.12	kN	Max Punch ;	53
Lower Punch :	Average ;	0.09	kN		
	Max ;	0.21	kN	Max Punch ;	47
Ejection Force :	Average ;	0.10	kN		
	Max ;	0.16	kN	Max Punch ;	22

Tablet Data (N/A)

Total Counter :	25792	Tablets			
A.W.C Counter :	23134	Tablets			
HEP :	76 (0.3%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	23057 (99.7%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	METHO 500MG STHA1A00
User ID :	kv190014	Batch NO. :	HAM2501150
Print Interval :	50000 Tablets	Running Time :	00 Hour 16 Min 28 Sec
Total Counter :	25792	Total Running Time :	2081 Hour 55 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	94	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		750000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		80	Min	Supply Time ;	0.5	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	1.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.28	kN	No		
HEP :	25.0 %	11.80	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	9.63	kN	Yes 20	Times	
Ref :		9.44	kN			
LCP :	2.0 %	9.25	kN	Yes 20	Times	
LEP :	25.0 %	7.08	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	6.61	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	0.77	kN	Min :	0.75	kN	Punch NO.;	41
SD :	2.50	%	Max :	0.88	kN	Punch NO.;	47
Mean Main Pressure :	9.35	kN	Min :	8.51	kN	Punch NO.;	1
SD :	4.90	%	Max :	10.35	kN	Punch NO.;	19
Filling Depth :	Increase ;	859	Times	Decrease ;	845	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	14	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	6.50	mm	Pre Thickness ;	4.50	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.66	mm	Main Thickness ;	2.66	mm
Filling Depth :	8.24	mm	Current Cam :	Cam C	
Main Air pressure :	0.792	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	9	Min			
Lower Punch Remain Time (S2):	31	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.07	kN		
	Max ;	0.14	kN	Max Punch ;	8
Lower Punch :	Average ;	0.09	kN		
	Max ;	0.21	kN	Max Punch ;	47
Ejection Force :	Average ;	0.10	kN		
	Max ;	0.16	kN	Max Punch ;	22

Tablet Data (N/A)

Total Counter :	25792	Tablets			
A.W.C Counter :	23134	Tablets			
HEP :	76 (0.3%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	23057 (99.7%)	Tablets			