

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Gaba800mg STJD1B00    |
| User ID :        | mr99799       | Batch NO. :          | JDM2501280            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 06 Min 45 Sec |
| Total Counter :  | 10942         | Total Running Time : | 2097 Hour 51 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |        |                        |         |                    |              |           |
|------------------------------------|--------|------------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 165                    | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam D<br>12.5-<br>18mm |         | Cam Range :        | Min :        | 10.50 mm  |
| Target Quantity :                  |        | 120000                 | Tablets |                    | Max :        | 16.00 mm  |
| Air Pressure Low Limit :           |        | 0.000                  | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                        |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                     | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                     | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                      | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00                   | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00                   | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00                   | kN      |                    |              |           |
|                                    |        |                        |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 24.91                  | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 23.95                  | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 19.55                  | kN      | Yes 15             | Times        |           |
| Ref :                              |        | 19.16                  | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 18.78                  | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 % | 14.37                  | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 30.0 % | 13.41                  | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                        |         | Yes                |              |           |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None   |                        |         |                    |              |           |

### Pressure Data

|                      |               |       |               |       |    |            |    |
|----------------------|---------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 2.20          | kN    | Min :         | 2.08  | kN | Punch NO.; | 45 |
| SD :                 | 2.50          | %     | Max :         | 2.30  | kN | Punch NO.; | 24 |
| Mean Main Pressure : | 19.35         | kN    | Min :         | 18.69 | kN | Punch NO.; | 37 |
| SD :                 | 1.80          | %     | Max :         | 20.06 | kN | Punch NO.; | 10 |
| Filling Depth :      | Increase ; 21 | Times | Decrease ; 23 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 16      | RPM | Capacity :       | 101760                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 26                     | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.01                   | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.60    | mm  | Main Thickness ; | 4.61                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 11.27   | mm  | Current Cam :    | Cam D<br>12.5-<br>18mm |              |
| Main Air pressure :           | 0.792   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 37      | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.27 | kN |             |    |
|                  | Max ;     | 0.30 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 3  |

### Tablet Data (N/A)

|                 |                  |         |       |        |         |
|-----------------|------------------|---------|-------|--------|---------|
| Total Counter : | 10942            | Tablets |       |        |         |
| A.W.C Counter : | 7897             | Tablets |       |        |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 7897<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Gaba800mg STJD1B00    |
| User ID :        | mr99799       | Batch NO. :          | JDM2501280            |
| Print Interval : | 50000 Tablets | Running Time :       | 00 Hour 06 Min 46 Sec |
| Total Counter :  | 10968         | Total Running Time : | 2097 Hour 51 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |        |                        |         |                    |              |           |
|------------------------------------|--------|------------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 180                    | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam D<br>12.5-<br>18mm |         | Cam Range :        | Min :        | 10.50 mm  |
| Target Quantity :                  |        | 120000                 | Tablets |                    | Max :        | 16.00 mm  |
| Air Pressure Low Limit :           |        | 0.000                  | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                        |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 80                     | Min     | Supply Time ;      | 0.5          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                     | Min     | Supply Time ;      | 1.0          | Sec       |
| Powder Supply Time :               |        | 60.0                   | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                      | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00                   | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00                   | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00                   | kN      |                    |              |           |
|                                    |        |                        |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 % | 25.56                  | kN      | No                 |              |           |
| HEP :                              | 25.0 % | 24.58                  | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 20.05                  | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 19.66                  | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 19.27                  | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 % | 14.75                  | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 % | 13.76                  | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                        |         | Yes                |              |           |
| Pre HSP                            |        | 40.00                  | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations              |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                        |         |                    |              |           |

### Pressure Data

|                      |              |       |               |       |    |            |    |
|----------------------|--------------|-------|---------------|-------|----|------------|----|
| Mean Pre Pressure :  | 1.98         | kN    | Min :         | 1.84  | kN | Punch NO.; | 16 |
| SD :                 | 2.90         | %     | Max :         | 2.09  | kN | Punch NO.; | 26 |
| Mean Main Pressure : | 19.90        | kN    | Min :         | 18.85 | kN | Punch NO.; | 28 |
| SD :                 | 1.80         | %     | Max :         | 20.75 | kN | Punch NO.; | 9  |
| Filling Depth :      | Increase ; 6 | Times | Decrease ; 12 | Times |    |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                        |              |
|-------------------------------|---------|-----|------------------|------------------------|--------------|
| Disk Speed :                  | 16      | RPM | Capacity :       | 101760                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 29                     | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |                        |              |
| Pre Lower Punch Position ;    | 8.00    | mm  | Pre Thickness ;  | 5.01                   | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |                        |              |
| Main Lower Punch Position ;   | 7.73    | mm  | Main Thickness ; | 4.74                   | mm           |
|                               |         |     |                  |                        |              |
| Filling Depth :               | 12.41   | mm  | Current Cam :    | Cam D<br>12.5-<br>18mm |              |
| Main Air pressure :           | 0.792   | MPa |                  |                        |              |
| Oil Lubrication status :      |         |     |                  |                        |              |
| Upper Punch Remain Time (S1): | 37      | Min |                  |                        |              |
| Lower Punch Remain Time (S2): | 19      | Min |                  |                        |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                        |              |
| Powder Status :               | Disable |     | Dust Collector : | Off                    |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off                    |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.27 | kN |             |    |
|                  | Max ;     | 0.30 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.11 | kN |             |    |
|                  | Max ;     | 0.19 | kN | Max Punch ; | 49 |
| Ejection Force : | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 3  |

### Tablet Data (N/A)

|                 |                  |         |       |        |         |
|-----------------|------------------|---------|-------|--------|---------|
| Total Counter : | 10968            | Tablets |       |        |         |
| A.W.C Counter : | 7897             | Tablets |       |        |         |
| HEP :           | 0 (0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 7897<br>(100.0%) | Tablets |       |        |         |