

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Tams0.4mg STUL2A00    |
| User ID :        | sj97904       | Batch NO. :          | ULM2501904            |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 02 Min 10 Sec |
| Total Counter :  | 1750786       | Total Running Time : | 2109 Hour 58 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 78       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2400000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 22.70    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 21.83    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 17.81    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 17.46    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 17.11    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 13.10    | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 12.22    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 1.71  | kN    | Min :      | 1.68  | kN    | Punch NO.; | 22 |
|                      |            | 0.70  | %     | Max :      | 1.74  | kN    | Punch NO.; | 28 |
| Mean Main Pressure : |            | 17.42 | kN    | Min :      | 15.24 | kN    | Punch NO.; | 19 |
|                      | SD :       | 4.60  | %     | Max :      | 18.83 | kN    | Punch NO.; | 49 |
| Filling Depth :      | Increase ; | 133   | Times | Decrease ; | 216   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.50    | mm  | Pre Thickness ;  | 5.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.29    | mm  | Main Thickness ; | 2.30     | mm           |
| Filling Depth :               | 8.19    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.805   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 6       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.29 | kN | Max Punch ; | 31 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 10 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 10 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1750786             | Tablets |       |        |         |
| A.W.C Counter : | 1741234             | Tablets |       |        |         |
| HEP :           | 13 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1741221<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | Tams0.4mg STUL2A00    |
| User ID :        | sj97904       | Batch NO. :          | ULM2501904            |
| Print Interval : | 50000 Tablets | Running Time :       | 11 Hour 02 Min 10 Sec |
| Total Counter :  | 1750786       | Total Running Time : | 2109 Hour 58 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 80       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2400000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 15       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 22.07    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 21.22    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 17.32    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 16.98    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 16.64    | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 12.73    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 11.88    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |                |       |            |       |       |            |    |
|----------------------|----------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 1.89           | kN    | Min :      | 1.83  | kN    | Punch NO.; | 1  |
| SD :                 | 1.20           | %     | Max :      | 1.96  | kN    | Punch NO.; | 8  |
| Mean Main Pressure : | 16.99          | kN    | Min :      | 15.07 | kN    | Punch NO.; | 19 |
| SD :                 | 5.10           | %     | Max :      | 19.08 | kN    | Punch NO.; | 45 |
| Filling Depth :      | Increase ; 147 | Times | Decrease ; | 188   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 25      | RPM | Capacity :       | 159000   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 20       | RPM          |
| Pre Upper Punch Position ;    | 2.99    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.50    | mm  | Pre Thickness ;  | 5.51     | mm           |
| Main Upper Punch Position ;   | 2.99    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.51    | mm  | Main Thickness ; | 2.52     | mm           |
| Filling Depth :               | 8.23    | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.805   | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 10      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 6       | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.23 | kN |             |    |
|                  | Max ;     | 0.29 | kN | Max Punch ; | 31 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.14 | kN | Max Punch ; | 10 |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 44 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1750786             | Tablets |       |        |         |
| A.W.C Counter : | 1741234             | Tablets |       |        |         |
| HEP :           | 13 (0.0%)           | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1741221<br>(100.0%) | Tablets |       |        |         |