

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 800MG STIU2D00
User ID :	sj97904	Batch NO. :	IUM2500980
Print Interval :	50000 Tablets	Running Time :	00 Hour 10 Min 58 Sec
Total Counter :	17156	Total Running Time :	2123 Hour 39 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	165	%	Manual ;	30	RPM
Filling Cam :		Cam D 12.5- 18mm		Cam Range :	Min :	10.50 mm
Target Quantity :		150000	Tablets		Max :	16.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	27.48	kN	No		
HEP :	25.0 %	26.42	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	21.56	kN	Yes 15	Times	
Ref :		21.14	kN			
LCP :	2.0 %	20.71	kN	Yes 15	Times	
LEP :	25.0 %	15.85	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	14.79	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.40	kN	Min :	3.30	kN	Punch NO.;	11
SD :	1.20	%	Max :	3.48	kN	Punch NO.;	47
Mean Main Pressure :	21.05	kN	Min :	19.44	kN	Punch NO.;	43
SD :	2.60	%	Max :	22.20	kN	Punch NO.;	51
Filling Depth :	Increase ; 66	Times	Decrease ; 48	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	25	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.70	mm	Main Thickness ;	4.71	mm
Filling Depth :	11.60	mm	Current Cam :	Cam D 12.5- 18mm	
Main Air pressure :	0.782	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.47	kN		
	Max ;	0.71	kN	Max Punch ;	31
Lower Punch :	Average ;	0.13	kN		
	Max ;	0.17	kN	Max Punch ;	50
Ejection Force :	Average ;	1.13	kN		
	Max ;	1.36	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	17156	Tablets			
A.W.C Counter :	15700	Tablets			
HEP :	0 (0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	15699 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 800MG STIU2D00
User ID :	sj97904	Batch NO. :	IUM2500980
Print Interval :	50000 Tablets	Running Time :	00 Hour 10 Min 58 Sec
Total Counter :	17156	Total Running Time :	2123 Hour 39 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	168	%	Manual ;	30	RPM
Filling Cam :		Cam D 12.5- 18mm		Cam Range :	Min :	10.50 mm
Target Quantity :		150000	Tablets		Max :	16.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	28.36	kN	No		
HEP :	25.0 %	27.26	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	22.25	kN	Yes 20	Times	
Ref :		21.81	kN			
LCP :	2.0 %	21.38	kN	Yes 20	Times	
LEP :	25.0 %	16.36	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	15.27	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.52	kN	Min :	3.38	kN	Punch NO.;	46
SD :	1.80	%	Max :	3.70	kN	Punch NO.;	30
Mean Main Pressure :	21.27	kN	Min :	19.73	kN	Punch NO.;	35
SD :	4.00	%	Max :	23.35	kN	Punch NO.;	51
Filling Depth :	Increase ;	112	Times	Decrease ;	118	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	25	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.85	mm	Main Thickness ;	4.86	mm
Filling Depth :	11.39	mm	Current Cam :	Cam D 12.5- 18mm	
Main Air pressure :	0.782	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	14	Min			
Lower Punch Remain Time (S2):	10	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.47	kN		
	Max ;	0.70	kN	Max Punch ;	31
Lower Punch :	Average ;	0.13	kN		
	Max ;	0.17	kN	Max Punch ;	50
Ejection Force :	Average ;	1.13	kN		
	Max ;	1.36	kN	Max Punch ;	17

Tablet Data (N/A)

Total Counter :	17156	Tablets			
A.W.C Counter :	15700	Tablets			
HEP :	0 (0%)	Tablets	LEP :	1 (0.0%)	Tablets
GOOD :	15699 (100.0%)	Tablets			