

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 800MG STIU2D00
User ID :	rp98927	Batch NO. :	IUM2501080
Print Interval :	50000 Tablets	Running Time :	00 Hour 03 Min 15 Sec
Total Counter :	5094	Total Running Time :	2125 Hour 4 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	170	%	Manual ;	30	RPM
Filling Cam :		Cam D 12.5- 18mm		Cam Range :	Min :	10.50 mm
Target Quantity :		150000	Tablets		Max :	16.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	28.26	kN	No		
HEP :	25.0 %	27.18	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	22.17	kN	Yes 15	Times	
Ref :		21.74	kN			
LCP :	2.0 %	21.31	kN	Yes 15	Times	
LEP :	25.0 %	16.31	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	15.22	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	3.44	kN	Min :	3.31	kN	Punch NO.;	14
SD :	1.50	%	Max :	3.60	kN	Punch NO.;	47
Mean Main Pressure :	20.86	kN	Min :	19.60	kN	Punch NO.;	32
SD :	2.60	%	Max :	21.84	kN	Punch NO.;	7
Filling Depth :	Increase ; 97	Times	Decrease ; 90	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	26	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.70	mm	Main Thickness ;	4.71	mm
Filling Depth :	11.74	mm	Current Cam :	Cam D 12.5- 18mm	
Main Air pressure :	0.790	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	12	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.41	kN		
	Max ;	0.63	kN	Max Punch ;	1
Lower Punch :	Average ;	0.15	kN		
	Max ;	0.19	kN	Max Punch ;	19
Ejection Force :	Average ;	1.32	kN		
	Max ;	1.51	kN	Max Punch ;	23

Tablet Data (N/A)

Total Counter :	5094	Tablets			
A.W.C Counter :	2738	Tablets			
HEP :	0 (0%)	Tablets	LEP :	23 (0.8%)	Tablets
GOOD :	2715 (99.2%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 800MG STIU2D00
User ID :	rp98927	Batch NO. :	IUM2501080
Print Interval :	50000 Tablets	Running Time :	00 Hour 03 Min 16 Sec
Total Counter :	5120	Total Running Time :	2125 Hour 4 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	182	%	Manual ;	30	RPM
Filling Cam :		Cam D 12.5- 18mm		Cam Range :	Min :	10.50 mm
Target Quantity :		150000	Tablets		Max :	16.00 mm
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		15	Min	Supply Time ;	2.0	Sec
Lower Punch Interval Time(S2) ;		15	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	29.72	kN	No		
HEP :	25.0 %	28.58	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	23.32	kN	Yes 20	Times	
Ref :		22.86	kN			
LCP :	2.0 %	22.41	kN	Yes 20	Times	
LEP :	25.0 %	17.15	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	16.00	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2	Rotations		Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	3.75	kN	Min :	3.62	kN	Punch NO.;	14
SD :	1.70	%	Max :	3.90	kN	Punch NO.;	42
Mean Main Pressure :	22.30	kN	Min :	20.05	kN	Punch NO.;	41
SD :	4.30	%	Max :	24.30	kN	Punch NO.;	2
Filling Depth :	Increase ;	200	Times	Decrease ;	213	Times	

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	95400	Tablets/hour
Feeder : Status ;	Auto		Speed ;	27	RPM
Pre Upper Punch Position ;	2.99	mm			
Pre Lower Punch Position ;	9.00	mm	Pre Thickness ;	6.01	mm
Main Upper Punch Position ;	2.99	mm			
Main Lower Punch Position ;	7.85	mm	Main Thickness ;	4.86	mm
Filling Depth :	11.50	mm	Current Cam :	Cam D 12.5- 18mm	
Main Air pressure :	0.789	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	12	Min			
Lower Punch Remain Time (S2):	8	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Enable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.43	kN		
	Max ;	0.65	kN	Max Punch ;	1
Lower Punch :	Average ;	0.15	kN		
	Max ;	0.19	kN	Max Punch ;	19
Ejection Force :	Average ;	1.32	kN		
	Max ;	1.51	kN	Max Punch ;	23

Tablet Data (N/A)

Total Counter :	5120	Tablets			
A.W.C Counter :	2738	Tablets			
HEP :	0 (0%)	Tablets	LEP :	23 (0.8%)	Tablets
GOOD :	2715 (99.2%)	Tablets			