

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 1000mg STYQ1A00   |
| User ID :        | pb97604       | Batch NO. :          | YQM2502110            |
| Print Interval : | 50000 Tablets | Running Time :       | 18 Hour 23 Min 26 Sec |
| Total Counter :  | 1850514       | Total Running Time : | 2165 Hour 16 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 147      | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 200      | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 90       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 26.48    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 25.47    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 20.78    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 20.37    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 19.96    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 15.28    | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 14.26    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 3.31  | kN    | Min :      | 3.12  | kN    | Punch NO.; | 17 |
|                      |            | 1.70  | %     | Max :      | 3.41  | kN    | Punch NO.; | 12 |
| Mean Main Pressure : |            | 20.18 | kN    | Min :      | 18.75 | kN    | Punch NO.; | 17 |
|                      | SD :       | 2.40  | %     | Max :      | 20.99 | kN    | Punch NO.; | 12 |
| Filling Depth :      | Increase ; | 448   | Times | Decrease ; | 437   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 16     | RPM | Capacity :       | 101760   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 24       | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 8.29   | mm  | Pre Thickness ;  | 5.30     | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 7.38   | mm  | Main Thickness ; | 4.39     | mm           |
| Filling Depth :               | 10.47  | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.808  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 44     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 36     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.07 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 21 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 1  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 1850514            | Tablets |       |             |         |
| A.W.C Counter : | 1847593            | Tablets |       |             |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 2619 (0.1%) | Tablets |
| GOOD :          | 1844974<br>(99.9%) | Tablets |       |             |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 1000mg STYQ1A00   |
| User ID :        | pb97604       | Batch NO. :          | YQM2502110            |
| Print Interval : | 50000 Tablets | Running Time :       | 18 Hour 23 Min 26 Sec |
| Total Counter :  | 1850514       | Total Running Time : | 2165 Hour 16 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 148      | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 200      | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 90       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 26.36    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 25.34    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 20.68    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 20.27    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 19.87    | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 15.21    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 14.19    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.44       | kN  | Min : | 3.26       | kN  | Punch NO.; | 35 |
| SD :                 | 2.30       | %   | Max : | 3.59       | kN  | Punch NO.; | 46 |
| Mean Main Pressure : | 19.90      | kN  | Min : | 18.38      | kN  | Punch NO.; | 5  |
| SD :                 | 3.30       | %   | Max : | 20.98      | kN  | Punch NO.; | 52 |
| Filling Depth :      | Increase ; | 272 | Times | Decrease ; | 270 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 16     | RPM | Capacity :       | 101760            | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 24                | RPM          |
| Pre Upper Punch Position ;    | 2.99   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 8.29   | mm  | Pre Thickness ;  | 5.30              | mm           |
| Main Upper Punch Position ;   | 2.99   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 7.55   | mm  | Main Thickness ; | 4.56              | mm           |
| Filling Depth :               | 10.26  | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.808  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 44     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 36     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.01 | kN |             |    |
|                  | Max ;     | 0.08 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.11 | kN | Max Punch ; | 21 |
| Ejection Force : | Average ; | 0.05 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 1  |

### Tablet Data (N/A)

|                 |                    |         |       |             |         |
|-----------------|--------------------|---------|-------|-------------|---------|
| Total Counter : | 1850514            | Tablets |       |             |         |
| A.W.C Counter : | 1847593            | Tablets |       |             |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 2619 (0.1%) | Tablets |
| GOOD :          | 1844974<br>(99.9%) | Tablets |       |             |         |