

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 1000mg STYQ1A00   |
| User ID :        | pb97604       | Batch NO. :          | YQM2502210            |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 19 Min 40 Sec |
| Total Counter :  | 1250364       | Total Running Time : | 2181 Hour 33 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 147      | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 25       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 29.05    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 27.93    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 22.79    | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 22.34    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 21.90    | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 16.76    | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 15.64    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |            |       |       |            |       |       |            |    |
|----------------------|------------|-------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | SD :       | 2.87  | kN    | Min :      | 2.77  | kN    | Punch NO.; | 35 |
|                      |            | 1.60  | %     | Max :      | 3.00  | kN    | Punch NO.; | 1  |
| Mean Main Pressure : | SD :       | 22.28 | kN    | Min :      | 21.43 | kN    | Punch NO.; | 19 |
|                      |            | 2.10  | %     | Max :      | 23.56 | kN    | Punch NO.; | 1  |
| Filling Depth :      | Increase ; | 547   | Times | Decrease ; | 563   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |                   |              |
|-------------------------------|--------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 16     | RPM | Capacity :       | 101760            | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 24                | RPM          |
| Pre Upper Punch Position ;    | 1.50   | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 7.00   | mm  | Pre Thickness ;  | 5.50              | mm           |
| Main Upper Punch Position ;   | 1.50   | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 6.15   | mm  | Main Thickness ; | 4.65              | mm           |
| Filling Depth :               | 10.10  | mm  | Current Cam :    | Cam C<br>8.5-14mm |              |
| Main Air pressure :           | 0.795  | MPa |                  |                   |              |
| Oil Lubrication status :      |        |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 64     | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 15     | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |                   |              |
| Powder Status :               | Enable |     | Dust Collector : | Off               |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.09 | kN | Max Punch ; | 27 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1250364             | Tablets |       |        |         |
| A.W.C Counter : | 1240972             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1240972<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MET 1000mg STYQ1A00   |
| User ID :        | pb97604       | Batch NO. :          | YQM2502210            |
| Print Interval : | 50000 Tablets | Running Time :       | 12 Hour 19 Min 40 Sec |
| Total Counter :  | 1250364       | Total Running Time : | 2181 Hour 33 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 148      | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam C    |         | Cam Range :        | Min :        | 6.50 mm   |
|                                    |             | 8.5-14mm |         |                    | Max :        | 12.00 mm  |
| Target Quantity :                  |             | 2250000  | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.400    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 25       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 29.09    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 27.97    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 22.83    | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 22.38    | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 21.93    | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 16.78    | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 15.66    | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |            |     |       |            |     |            |    |
|----------------------|------------|-----|-------|------------|-----|------------|----|
| Mean Pre Pressure :  | 3.06       | kN  | Min : | 2.92       | kN  | Punch NO.; | 35 |
| SD :                 | 2.00       | %   | Max : | 3.23       | kN  | Punch NO.; | 12 |
| Mean Main Pressure : | 22.20      | kN  | Min : | 21.05      | kN  | Punch NO.; | 35 |
| SD :                 | 2.90       | %   | Max : | 23.72      | kN  | Punch NO.; | 12 |
| Filling Depth :      | Increase ; | 348 | Times | Decrease ; | 346 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |        |     |                  |          |              |
|-------------------------------|--------|-----|------------------|----------|--------------|
| Disk Speed :                  | 16     | RPM | Capacity :       | 101760   | Tablets/hour |
| Feeder : Status ;             | Auto   |     | Speed ;          | 24       | RPM          |
| Pre Upper Punch Position ;    | 1.50   | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 7.00   | mm  | Pre Thickness ;  | 5.50     | mm           |
| Main Upper Punch Position ;   | 1.50   | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.20   | mm  | Main Thickness ; | 4.70     | mm           |
| Filling Depth :               | 9.98   | mm  | Current Cam :    | Cam C    |              |
| Main Air pressure :           | 0.795  | MPa |                  | 8.5-14mm |              |
| Oil Lubrication status :      |        |     |                  |          |              |
| Upper Punch Remain Time (S1): | 64     | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 15     | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0      | Min |                  |          |              |
| Powder Status :               | Enable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off    |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.02 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 25 |
| Lower Punch :    | Average ; | 0.09 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 1  |
| Ejection Force : | Average ; | 0.04 | kN |             |    |
|                  | Max ;     | 0.10 | kN | Max Punch ; | 46 |

### Tablet Data (N/A)

|                 |                     |         |       |        |         |
|-----------------|---------------------|---------|-------|--------|---------|
| Total Counter : | 1250364             | Tablets |       |        |         |
| A.W.C Counter : | 1240972             | Tablets |       |        |         |
| HEP :           | 0 (0%)              | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 1240972<br>(100.0%) | Tablets |       |        |         |