

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 400MG STIU2B00
User ID :	pb97604	Batch NO. :	IUM2500740
Print Interval :	50000 Tablets	Running Time :	03 Hour 35 Min 24 Sec
Total Counter :	454874	Total Running Time :	2194 Hour 59 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	241	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	27.84	kN	No		
HEP :	25.0 %	26.77	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	21.85	kN	Yes 15	Times	
Ref :		21.42	kN			
LCP :	2.0 %	20.99	kN	Yes 15	Times	
LEP :	25.0 %	16.06	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	14.99	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.52	kN	Min :	2.41	kN	Punch NO.;	52
SD :	3.10	%	Max :	2.73	kN	Punch NO.;	21
Mean Main Pressure :	20.12	kN	Min :	18.29	kN	Punch NO.;	52
SD :	6.70	%	Max :	23.72	kN	Punch NO.;	21
Filling Depth :	Increase ; 121	Times	Decrease ;	78	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	48	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	5.20	mm	Pre Thickness ;	3.20	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.40	mm	Main Thickness ;	2.40	mm
Filling Depth :	9.79	mm	Current Cam :	Cam C 8.5-14mm	
Main Air pressure :	0.800	MPa			
Oil Lubrication status :					
Upper Punch Remain Time (S1):	43	Min			
Lower Punch Remain Time (S2):	11	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.55	kN	Max Punch ;	12
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.07	kN	Max Punch ;	24

Tablet Data (N/A)

Total Counter :	454874	Tablets			
A.W.C Counter :	448504	Tablets			
HEP :	626 (0.1%)	Tablets	LEP :	9 (0.0%)	Tablets
GOOD :	447869 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	IBU 400MG STIU2B00
User ID :	pb97604	Batch NO. :	IUM2500740
Print Interval :	50000 Tablets	Running Time :	03 Hour 35 Min 24 Sec
Total Counter :	454874	Total Running Time :	2194 Hour 59 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	242	%	Manual ;	30	RPM
Filling Cam :		Cam C		Cam Range :	Min :	6.50 mm
		8.5-14mm			Max :	12.00 mm
Target Quantity :		2025000	Tablets			
Air Pressure Low Limit :		0.400	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	30.43	kN	No		
HEP :	25.0 %	29.26	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	23.88	kN	Yes 20	Times	
Ref :		23.41	kN			
LCP :	2.0 %	22.94	kN	Yes 20	Times	
LEP :	25.0 %	17.56	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	16.39	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	SD :	2.76	kN	Min :	2.63	kN	Punch NO.;	52
		3.00	%	Max :	2.93	kN	Punch NO.;	25
Mean Main Pressure :		23.07	kN	Min :	20.59	kN	Punch NO.;	31
	SD :	6.60	%	Max :	26.79	kN	Punch NO.;	48
Filling Depth :	Increase ;	88	Times	Decrease ;	73	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	20	RPM	Capacity :	0	Tablets/hour
Feeder : Status ;	Auto		Speed ;	48	RPM
Pre Upper Punch Position ;	2.00	mm			
Pre Lower Punch Position ;	5.20	mm	Pre Thickness ;	3.20	mm
Main Upper Punch Position ;	2.00	mm			
Main Lower Punch Position ;	4.40	mm	Main Thickness ;	2.40	mm
Filling Depth :	8.86	mm	Current Cam :	Cam C	
Main Air pressure :	0.800	MPa		8.5-14mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	43	Min			
Lower Punch Remain Time (S2):	11	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	On	

Tightness

Upper Punch :	Average ;	0.00	kN		
	Max ;	0.00	kN	Max Punch ;	0
Lower Punch :	Average ;	0.14	kN		
	Max ;	0.51	kN	Max Punch ;	12
Ejection Force :	Average ;	0.04	kN		
	Max ;	0.07	kN	Max Punch ;	24

Tablet Data (N/A)

Total Counter :	454874	Tablets			
A.W.C Counter :	448504	Tablets			
HEP :	626 (0.1%)	Tablets	LEP :	9 (0.0%)	Tablets
GOOD :	447869 (99.9%)	Tablets			