

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MIRTAZAPINE 15 MG
User ID :	sb98304	Batch NO. :	MRT TAB XV 25 001B
Print Interval :	50000 Tablets	Running Time :	00 Hour 21 Min 59 Sec
Total Counter :	50112	Total Running Time :	2207 Hour 36 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	20	%	Manual ;	57	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		513333	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.20	kN	No		
HEP :	25.0 %	11.73	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	9.57	kN	Yes 15	Times	
Ref :		9.38	kN			
LCP :	2.0 %	9.19	kN	Yes 15	Times	
LEP :	25.0 %	7.04	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	6.57	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.36	kN	Min :	2.30	kN	Punch NO.;	47
SD :	1.00	%	Max :	2.42	kN	Punch NO.;	34
Mean Main Pressure :	9.07	kN	Min :	8.86	kN	Punch NO.;	47
SD :	0.90	%	Max :	9.27	kN	Punch NO.;	57
Filling Depth :	Increase ; 1	Times	Decrease ;	9	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	231000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	5	RPM
Pre Upper Punch Position ;	4.50	mm			
Pre Lower Punch Position ;	6.65	mm	Pre Thickness ;	2.15	mm
Main Upper Punch Position ;	4.50	mm			
Main Lower Punch Position ;	6.35	mm	Main Thickness ;	1.85	mm
Filling Depth :	7.70	mm	Current Cam :	Cam B	
Main Air pressure :	0.792	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	2	Min			
Lower Punch Remain Time (S2):	30	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.18	kN		
	Max ;	0.21	kN	Max Punch ;	53
Lower Punch :	Average ;	0.07	kN		
	Max ;	0.13	kN	Max Punch ;	7
Ejection Force :	Average ;	0.16	kN		
	Max ;	0.21	kN	Max Punch ;	46

Tablet Data (N/A)

Total Counter :	50112	Tablets			
A.W.C Counter :	35376	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	35376 (100.0%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name : VANTIX P720D User ID : sb98304 Print Interval : 50000 Tablets Total Counter : 50176 Station : Station 2	Product Name : MIRTAZAPINE 15 MG Batch NO. : MRT TAB XV 25 001B Running Time : 00 Hour 22 Min 00 Sec Total Running Time : 2207 Hour 36 Min
---	---

Setting Value

Feeder :	Auto ;	20	%	Manual ;	30	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
Target Quantity :		4.5-10mm			Max :	8.00 mm
Air Pressure Low Limit :		513333	Tablets			
Oil Lubrication :		0.000	MPa			
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	12.00	kN	No		
HEP :	25.0 %	11.54	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	9.41	kN	Yes 20	Times	
Ref :		9.23	kN			
LCP :	2.0 %	9.05	kN	Yes 20	Times	
LEP :	25.0 %	6.92	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	6.46	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.53	kN	Min :	2.46	kN	Punch NO.;	4
SD :	1.40	%	Max :	2.61	kN	Punch NO.;	69
Mean Main Pressure :	9.70	kN	Min :	9.12	kN	Punch NO.;	27
SD :	2.80	%	Max :	10.38	kN	Punch NO.;	72
Filling Depth :	Increase ; 1	Times	Decrease ; 1	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	25	RPM	Capacity :	231000	Tablets/hour
Feeder : Status ;	Auto		Speed ;	5	RPM
Pre Upper Punch Position ;	4.51	mm			
Pre Lower Punch Position ;	6.80	mm	Pre Thickness ;	2.29	mm
Main Upper Punch Position ;	4.50	mm			
Main Lower Punch Position ;	6.55	mm	Main Thickness ;	2.05	mm
Filling Depth :	6.70	mm	Current Cam :	Cam B	
Main Air pressure :	0.790	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	2	Min			
Lower Punch Remain Time (S2):	30	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	On		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.18	kN		
	Max ;	0.21	kN	Max Punch ;	53
Lower Punch :	Average ;	0.07	kN		
	Max ;	0.13	kN	Max Punch ;	54
Ejection Force :	Average ;	0.16	kN		
	Max ;	0.21	kN	Max Punch ;	49

Tablet Data (N/A)

Total Counter :	50176	Tablets			
A.W.C Counter :	35376	Tablets			
HEP :	0 (0%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	35376 (100.0%)	Tablets			