

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MIRTAZAPINE 15 MG     |
| User ID :        | rp98927       | Batch NO. :          | MRT TAB XV 25 001B    |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 36 Min 49 Sec |
| Total Counter :  | 300166        | Total Running Time : | 2208 Hour 50 Min      |
| Station :        | Station 1     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 20       | %       | Manual ;           | 57           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 513333   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 12.27    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 11.80    | kN      | Yes 10             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 9.63     | kN      | Yes 15             | Times        |           |
| Ref :                              |             | 9.44     | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 9.25     | kN      | Yes 15             | Times        |           |
| LEP :                              | 25.0 %      | 7.08     | kN      | No 15              | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 6.61     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | None        |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |      |       |            |    |
|----------------------|--------------|-------|------------|------|-------|------------|----|
| Mean Pre Pressure :  | 2.44         | kN    | Min :      | 2.36 | kN    | Punch NO.; | 12 |
| SD :                 | 1.00         | %     | Max :      | 2.50 | kN    | Punch NO.; | 2  |
| Mean Main Pressure : | 9.48         | kN    | Min :      | 9.11 | kN    | Punch NO.; | 12 |
| SD :                 | 1.50         | %     | Max :      | 9.78 | kN    | Punch NO.; | 45 |
| Filling Depth :      | Increase ; 5 | Times | Decrease ; | 12   | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0        | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4        | RPM          |
| Pre Upper Punch Position ;    | 4.50    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.65    | mm  | Pre Thickness ;  | 2.15     | mm           |
| Main Upper Punch Position ;   | 4.50    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 6.35    | mm  | Main Thickness ; | 1.85     | mm           |
| Filling Depth :               | 6.63    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.770   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 40      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | On      |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.20 | kN | Max Punch ; | 72 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 56 |
| Ejection Force : | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 27 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 300166            | Tablets |       |        |         |
| A.W.C Counter : | 276090            | Tablets |       |        |         |
| HEP :           | 376 (0.1%)        | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 275714<br>(99.9%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                       |
|------------------|---------------|----------------------|-----------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | MIRTAZAPINE 15 MG     |
| User ID :        | rp98927       | Batch NO. :          | MRT TAB XV 25 001B    |
| Print Interval : | 50000 Tablets | Running Time :       | 01 Hour 36 Min 49 Sec |
| Total Counter :  | 300166        | Total Running Time : | 2208 Hour 50 Min      |
| Station :        | Station 2     |                      |                       |

### Setting Value

|                                    |             |          |         |                    |              |           |
|------------------------------------|-------------|----------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ;      | 20       | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |             | Cam B    |         | Cam Range :        | Min :        | 2.50 mm   |
|                                    |             | 4.5-10mm |         |                    | Max :        | 8.00 mm   |
| Target Quantity :                  |             | 513333   | Tablets |                    |              |           |
| Air Pressure Low Limit :           |             | 0.000    | MPa     |                    |              |           |
| Oil Lubrication :                  |             |          |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |             | 70       | Min     | Supply Time ;      | 1.0          | Sec       |
| Lower Punch Interval Time(S2) ;    |             | 40       | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |             | 60.0     | Sec     |                    |              |           |
| Initial Reject Time :              |             | 2        | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |             | 1.00     | kN      |                    |              |           |
| Ejecting Force High Limit :        |             | 4.00     | kN      |                    |              |           |
|                                    |             |          |         | Stop Condition ;   |              |           |
| HSP :                              | 30.0 %      | 12.40    | kN      | No                 |              |           |
| HEP :                              | 25.0 %      | 11.93    | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %       | 9.73     | kN      | Yes 20             | Times        |           |
| Ref :                              |             | 9.54     | kN      |                    |              |           |
| LCP :                              | 2.0 %       | 9.35     | kN      | Yes 20             | Times        |           |
| LEP :                              | 25.0 %      | 7.16     | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 30.0 %      | 6.68     | kN      | No                 |              |           |
| SD limit :                         | 45.0 %      |          |         | Yes                |              |           |
| Pre HSP                            |             | 40.00    | kN      | Yes                |              |           |
| Mean Calculation range :           | 2 Rotations |          |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N           |          |         |                    |              |           |

### Pressure Data

|                      |              |       |            |       |       |            |    |
|----------------------|--------------|-------|------------|-------|-------|------------|----|
| Mean Pre Pressure :  | 2.47         | kN    | Min :      | 2.41  | kN    | Punch NO.; | 15 |
| SD :                 | 1.00         | %     | Max :      | 2.53  | kN    | Punch NO.; | 25 |
| Mean Main Pressure : | 9.51         | kN    | Min :      | 8.92  | kN    | Punch NO.; | 63 |
| SD :                 | 2.70         | %     | Max :      | 10.23 | kN    | Punch NO.; | 48 |
| Filling Depth :      | Increase ; 5 | Times | Decrease ; | 27    | Times |            |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |                   |              |
|-------------------------------|---------|-----|------------------|-------------------|--------------|
| Disk Speed :                  | 20      | RPM | Capacity :       | 0                 | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 4                 | RPM          |
| Pre Upper Punch Position ;    | 4.51    | mm  |                  |                   |              |
| Pre Lower Punch Position ;    | 6.80    | mm  | Pre Thickness ;  | 2.29              | mm           |
| Main Upper Punch Position ;   | 4.50    | mm  |                  |                   |              |
| Main Lower Punch Position ;   | 6.52    | mm  | Main Thickness ; | 2.02              | mm           |
| Filling Depth :               | 6.34    | mm  | Current Cam :    | Cam B<br>4.5-10mm |              |
| Main Air pressure :           | 0.770   | MPa |                  |                   |              |
| Oil Lubrication status :      |         |     |                  |                   |              |
| Upper Punch Remain Time (S1): | 2       | Min |                  |                   |              |
| Lower Punch Remain Time (S2): | 40      | Min |                  |                   |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |                   |              |
| Powder Status :               | Disable |     | Dust Collector : | Off               |              |
| Initial Reject :              | On      |     | Buzzer :         | Off               |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.17 | kN |             |    |
|                  | Max ;     | 0.21 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.10 | kN |             |    |
|                  | Max ;     | 0.13 | kN | Max Punch ; | 56 |
| Ejection Force : | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.16 | kN | Max Punch ; | 27 |

### Tablet Data (N/A)

|                 |                   |         |       |        |         |
|-----------------|-------------------|---------|-------|--------|---------|
| Total Counter : | 300166            | Tablets |       |        |         |
| A.W.C Counter : | 276090            | Tablets |       |        |         |
| HEP :           | 376 (0.1%)        | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 275714<br>(99.9%) | Tablets |       |        |         |