

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MIRTAZAPINE 15 MG
User ID :	sj97904	Batch NO. :	MRT TAB XV 25 001B
Print Interval :	50000 Tablets	Running Time :	01 Hour 58 Min 19 Sec
Total Counter :	350182	Total Running Time :	2209 Hour 12 Min
Station :	Station 1		

Setting Value

Feeder :	Auto ;	22	%	Manual ;	57	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		513333	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	13.53	kN	No		
HEP :	25.0 %	13.01	kN	Yes 10	Rotations 20	Tablets
HCP :	2.0 %	10.62	kN	Yes 15	Times	
Ref :		10.41	kN			
LCP :	2.0 %	10.20	kN	Yes 15	Times	
LEP :	25.0 %	7.81	kN	No 15	Rotations 20	Tablets
LSP :	30.0 %	7.29	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	None					

Pressure Data

Mean Pre Pressure :	2.84	kN	Min :	2.69	kN	Punch NO.;	18
SD :	2.20	%	Max :	3.02	kN	Punch NO.;	32
Mean Main Pressure :	10.75	kN	Min :	10.07	kN	Punch NO.;	18
SD :	3.30	%	Max :	11.73	kN	Punch NO.;	32
Filling Depth :	Increase ; 25	Times	Decrease ; 26	Times			

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	138600	Tablets/hour
Feeder : Status ;	Auto		Speed ;	3	RPM
Pre Upper Punch Position ;	4.50	mm			
Pre Lower Punch Position ;	6.30	mm	Pre Thickness ;	1.80	mm
Main Upper Punch Position ;	4.50	mm			
Main Lower Punch Position ;	6.20	mm	Main Thickness ;	1.70	mm
Filling Depth :	6.54	mm	Current Cam :	Cam B	
Main Air pressure :	0.676	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	55	Min			
Lower Punch Remain Time (S2):	23	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.15	kN		
	Max ;	0.17	kN	Max Punch ;	33
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.10	kN	Max Punch ;	58
Ejection Force :	Average ;	0.12	kN		
	Max ;	0.15	kN	Max Punch ;	31

Tablet Data (N/A)

Total Counter :	350182	Tablets			
A.W.C Counter :	318492	Tablets			
HEP :	376 (0.1%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	318116 (99.9%)	Tablets			

PRODUCTION REPORT

Product Information

Machine Name :	VANTIX P720D	Product Name :	MIRTAZAPINE 15 MG
User ID :	sj97904	Batch NO. :	MRT TAB XV 25 001B
Print Interval :	50000 Tablets	Running Time :	01 Hour 58 Min 19 Sec
Total Counter :	350182	Total Running Time :	2209 Hour 12 Min
Station :	Station 2		

Setting Value

Feeder :	Auto ;	20	%	Manual ;	30	RPM
Filling Cam :		Cam B		Cam Range :	Min :	2.50 mm
		4.5-10mm			Max :	8.00 mm
Target Quantity :		513333	Tablets			
Air Pressure Low Limit :		0.000	MPa			
Oil Lubrication :						
Upper Punch Interval Time(S1) ;		70	Min	Supply Time ;	1.0	Sec
Lower Punch Interval Time(S2) ;		40	Min	Supply Time ;	2.0	Sec
Powder Supply Time :		60.0	Sec			
Initial Reject Time :		2	Tablets			
Upper Punch Tightness High Limit :		1.00	kN			
Lower Punch Tightness High Limit :		1.00	kN			
Ejecting Force High Limit :		4.00	kN			
				Stop Condition ;		
HSP :	30.0 %	13.77	kN	No		
HEP :	25.0 %	13.24	kN	Yes 15	Rotations 20	Tablets
HCP :	2.0 %	10.81	kN	Yes 20	Times	
Ref :		10.59	kN			
LCP :	2.0 %	10.38	kN	Yes 20	Times	
LEP :	25.0 %	7.95	kN	Yes 15	Rotations 20	Tablets
LSP :	30.0 %	7.42	kN	No		
SD limit :	45.0 %			Yes		
Pre HSP		40.00	kN	Yes		
Mean Calculation range :	2 Rotations			Calculation Pass :	2	Rotations
Empty Punch No. :	N					

Pressure Data

Mean Pre Pressure :	2.86	kN	Min :	2.80	kN	Punch NO.;	3
SD :	0.90	%	Max :	2.92	kN	Punch NO.;	69
Mean Main Pressure :	10.61	kN	Min :	10.06	kN	Punch NO.;	63
SD :	2.00	%	Max :	11.18	kN	Punch NO.;	69
Filling Depth :	Increase ; 5	Times	Decrease ;	27	Times		

PRODUCTION REPORT

Operation Value

Disk Speed :	15	RPM	Capacity :	138600	Tablets/hour
Feeder : Status ;	Auto		Speed ;	3	RPM
Pre Upper Punch Position ;	4.51	mm			
Pre Lower Punch Position ;	6.36	mm	Pre Thickness ;	1.85	mm
Main Upper Punch Position ;	4.50	mm			
Main Lower Punch Position ;	6.40	mm	Main Thickness ;	1.90	mm
Filling Depth :	6.23	mm	Current Cam :	Cam B	
Main Air pressure :	0.676	MPa		4.5-10mm	
Oil Lubrication status :					
Upper Punch Remain Time (S1):	55	Min			
Lower Punch Remain Time (S2):	23	Min			
Lower Head Remain Time (S3):	0	Min			
Powder Status :	Disable		Dust Collector :	Off	
Initial Reject :	Off		Buzzer :	Off	

Tightness

Upper Punch :	Average ;	0.15	kN		
	Max ;	0.17	kN	Max Punch ;	33
Lower Punch :	Average ;	0.08	kN		
	Max ;	0.10	kN	Max Punch ;	58
Ejection Force :	Average ;	0.12	kN		
	Max ;	0.15	kN	Max Punch ;	19

Tablet Data (N/A)

Total Counter :	350182	Tablets			
A.W.C Counter :	318492	Tablets			
HEP :	376 (0.1%)	Tablets	LEP :	0 (0%)	Tablets
GOOD :	318116 (99.9%)	Tablets			