

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | LER HCL20MG<br>STDJ1B00 |
| User ID :        | rp98927       | Batch NO. :          | DJM2500520              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 57 Min 39 Sec   |
| Total Counter :  | 400106        | Total Running Time : | 2257 Hour 28 Min        |
| Station :        | Station 1     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |       |                      |
|------------------------------------|--------|-------------------|---------|--------------------|-------|----------------------|
| Feeder :                           | Auto ; | 50                | %       | Manual ;           | 57    | RPM                  |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min : | 2.50 mm              |
| Target Quantity :                  |        | 2550000           | Tablets |                    | Max : | 8.00 mm              |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |       |                      |
| Oil Lubrication :                  |        |                   |         |                    |       |                      |
| Upper Punch Interval Time(S1) ;    |        | 35                | Min     | Supply Time ;      | 0.3   | Sec                  |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 2.0   | Sec                  |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |       |                      |
| Initial Reject Time :              |        | 2                 | Tablets |                    |       |                      |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |       |                      |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |       |                      |
|                                    |        |                   |         | Stop Condition ;   |       |                      |
| HSP :                              | 25.0 % | 12.85             | kN      | No                 |       |                      |
| HEP :                              | 20.0 % | 12.33             | kN      | Yes                | 10    | Rotations 20 Tablets |
| HCP :                              | 2.0 %  | 10.48             | kN      | Yes                | 15    | Times                |
| Ref :                              |        | 10.28             | kN      |                    |       |                      |
| LCP :                              | 2.0 %  | 10.07             | kN      | Yes                | 15    | Times                |
| LEP :                              | 20.0 % | 8.22              | kN      | No                 | 15    | Rotations 20 Tablets |
| LSP :                              | 25.0 % | 7.71              | kN      | No                 |       |                      |
| SD limit :                         | 45.0 % |                   |         | Yes                |       |                      |
| Pre HSP                            |        | 40.00             | kN      | Yes                |       |                      |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2     | Rotations            |
| Empty Punch No. :                  | None   |                   |         |                    |       |                      |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.55       | kN | Min : | 2.50       | kN | Punch NO.; | 20 |
| SD :                 | 1.00       | %  | Max : | 2.61       | kN | Punch NO.; | 31 |
| Mean Main Pressure : | 10.33      | kN | Min : | 9.64       | kN | Punch NO.; | 41 |
| SD :                 | 2.80       | %  | Max : | 11.11      | kN | Punch NO.; | 27 |
| Filling Depth :      | Increase ; | 6  | Times | Decrease ; | 10 | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 11       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.05    | mm  | Main Thickness ; | 2.05     | mm           |
| Filling Depth :               | 5.50    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.822   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 34      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 17      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.33 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.17 | kN | Max Punch ; | 24 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 50 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 400106             | Tablets |       |        |         |
| A.W.C Counter : | 398557             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 398557<br>(100.0%) | Tablets |       |        |         |

## PRODUCTION REPORT

### Product Information

|                  |               |                      |                         |
|------------------|---------------|----------------------|-------------------------|
| Machine Name :   | VANTIX P720D  | Product Name :       | LER HCL20MG<br>STDJ1B00 |
| User ID :        | rp98927       | Batch NO. :          | DJM2500520              |
| Print Interval : | 50000 Tablets | Running Time :       | 02 Hour 57 Min 39 Sec   |
| Total Counter :  | 400106        | Total Running Time : | 2257 Hour 28 Min        |
| Station :        | Station 2     |                      |                         |

### Setting Value

|                                    |        |                   |         |                    |              |           |
|------------------------------------|--------|-------------------|---------|--------------------|--------------|-----------|
| Feeder :                           | Auto ; | 50                | %       | Manual ;           | 30           | RPM       |
| Filling Cam :                      |        | Cam B<br>4.5-10mm |         | Cam Range :        | Min :        | 2.50 mm   |
| Target Quantity :                  |        | 2550000           | Tablets |                    | Max :        | 8.00 mm   |
| Air Pressure Low Limit :           |        | 0.000             | MPa     |                    |              |           |
| Oil Lubrication :                  |        |                   |         |                    |              |           |
| Upper Punch Interval Time(S1) ;    |        | 35                | Min     | Supply Time ;      | 0.3          | Sec       |
| Lower Punch Interval Time(S2) ;    |        | 40                | Min     | Supply Time ;      | 2.0          | Sec       |
| Powder Supply Time :               |        | 60.0              | Sec     |                    |              |           |
| Initial Reject Time :              |        | 2                 | Tablets |                    |              |           |
| Upper Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Lower Punch Tightness High Limit : |        | 1.00              | kN      |                    |              |           |
| Ejecting Force High Limit :        |        | 4.00              | kN      |                    |              |           |
|                                    |        |                   |         | Stop Condition ;   |              |           |
| HSP :                              | 25.0 % | 12.84             | kN      | No                 |              |           |
| HEP :                              | 20.0 % | 12.33             | kN      | Yes 15             | Rotations 20 | Tablets   |
| HCP :                              | 2.0 %  | 10.48             | kN      | Yes 20             | Times        |           |
| Ref :                              |        | 10.27             | kN      |                    |              |           |
| LCP :                              | 2.0 %  | 10.07             | kN      | Yes 20             | Times        |           |
| LEP :                              | 20.0 % | 8.22              | kN      | Yes 15             | Rotations 20 | Tablets   |
| LSP :                              | 25.0 % | 7.70              | kN      | No                 |              |           |
| SD limit :                         | 45.0 % |                   |         | Yes                |              |           |
| Pre HSP                            |        | 40.00             | kN      | Yes                |              |           |
| Mean Calculation range :           | 2      | Rotations         |         | Calculation Pass : | 2            | Rotations |
| Empty Punch No. :                  | N      |                   |         |                    |              |           |

### Pressure Data

|                      |            |    |       |            |    |            |    |
|----------------------|------------|----|-------|------------|----|------------|----|
| Mean Pre Pressure :  | 2.75       | kN | Min : | 2.68       | kN | Punch NO.; | 13 |
| SD :                 | 1.30       | %  | Max : | 2.83       | kN | Punch NO.; | 14 |
| Mean Main Pressure : | 10.34      | kN | Min : | 9.71       | kN | Punch NO.; | 15 |
| SD :                 | 2.50       | %  | Max : | 11.15      | kN | Punch NO.; | 27 |
| Filling Depth :      | Increase ; | 4  | Times | Decrease ; | 2  | Times      |    |

## PRODUCTION REPORT

### Operation Value

|                               |         |     |                  |          |              |
|-------------------------------|---------|-----|------------------|----------|--------------|
| Disk Speed :                  | 22      | RPM | Capacity :       | 139920   | Tablets/hour |
| Feeder : Status ;             | Auto    |     | Speed ;          | 11       | RPM          |
| Pre Upper Punch Position ;    | 3.00    | mm  |                  |          |              |
| Pre Lower Punch Position ;    | 6.00    | mm  | Pre Thickness ;  | 3.00     | mm           |
| Main Upper Punch Position ;   | 3.00    | mm  |                  |          |              |
| Main Lower Punch Position ;   | 5.30    | mm  | Main Thickness ; | 2.30     | mm           |
| Filling Depth :               | 6.06    | mm  | Current Cam :    | Cam B    |              |
| Main Air pressure :           | 0.822   | MPa |                  | 4.5-10mm |              |
| Oil Lubrication status :      |         |     |                  |          |              |
| Upper Punch Remain Time (S1): | 34      | Min |                  |          |              |
| Lower Punch Remain Time (S2): | 17      | Min |                  |          |              |
| Lower Head Remain Time (S3):  | 0       | Min |                  |          |              |
| Powder Status :               | Disable |     | Dust Collector : | Off      |              |
| Initial Reject :              | Off     |     | Buzzer :         | Off      |              |

### Tightness

|                  |           |      |    |             |    |
|------------------|-----------|------|----|-------------|----|
| Upper Punch :    | Average ; | 0.26 | kN |             |    |
|                  | Max ;     | 0.33 | kN | Max Punch ; | 35 |
| Lower Punch :    | Average ; | 0.12 | kN |             |    |
|                  | Max ;     | 0.15 | kN | Max Punch ; | 15 |
| Ejection Force : | Average ; | 0.03 | kN |             |    |
|                  | Max ;     | 0.05 | kN | Max Punch ; | 50 |

### Tablet Data (N/A)

|                 |                    |         |       |        |         |
|-----------------|--------------------|---------|-------|--------|---------|
| Total Counter : | 400106             | Tablets |       |        |         |
| A.W.C Counter : | 398557             | Tablets |       |        |         |
| HEP :           | 0 (0%)             | Tablets | LEP : | 0 (0%) | Tablets |
| GOOD :          | 398557<br>(100.0%) | Tablets |       |        |         |